

Support Guide



Introduction

This Support Guide is intended to help clinicians understand and use the *GI Effects*® *Comprehensive Profile*, a select set of fecal biomarkers aimed at identifying key processes that influence both gastrointestinal and overall health.

Table of Contents



GI *fx* GI Effects Stool Profiles®

- Interpretation-At-A-Glance Synthesis3
 - Four Functional Pillars
 - Diversity and Abundance
- Biomarker Review.....4
 - Organization of the Biomarker Review
 - Graphical Representation and Color Coding for Biomarkers
- Digestion and Absorption7
 - Pancreatic Elastase-1 (PE1)
 - Products of Protein Breakdown
 - Total Fecal Fats
- Inflammation and Immunology11
 - Calprotectin
 - Eosinophil Protein X (EPX)
 - Fecal Secretory IgA (sIgA)
 - Lactoferrin (Add-on)
- Gastrointestinal Microbiome.....16
 - Metabolic Products: Short-Chain Fatty Acids
 - Beta-glucuronidase
 - Commensal Bacteria
 - Individual Commensal Bacteria
 - Firmicutes/Bacteroidetes Ratio
 - Bacteriology (Culture)
 - Additional Bacteria (Culture)
 - Mycology (Culture)
 - Parasitology

- Common Parasitic Protozoans30
- Bacteria Sensitivity.....31
- Mycology Sensitivity.....31
- Additional Tests.....31
 - Color
 - Consistency
 - Occult Blood
 - Pathogenic Bacteria (Add-on)
- Pathogenic Bacteria by EIA33
- References34
- Gut Microbiota Clinical Association Summary.....43

Leading the Evolution in Gut Health Assessment

Interpretation-At-A-Glance Synthesis

The test report is organized so that the clinician may move through the results in a logical order that enhances clinical utility, beginning with the Interpretation-At-A-Glance Synthesis Page.


Using evidence-based rules and weighted algorithms, this page synthesizes patient test results into key functional areas of clinical significance and provides a directional indication of potential next steps in patient management.

The Interpretation-At-A-Glance page is divided into two major sections: 1) Four Functional Pillars, and 2) Diversity and Abundance.

Four Functional Pillars

Pertinent biomarkers have been grouped into four clinically actionable areas:

- Infection
- Inflammation
- Insufficiency (Digestive)
- Imbalance (Metabolic)



The four functional pillars utilize a proprietary algorithm to evaluate key clinical markers in the four functional areas. The algorithm takes into account the level of each individual biomarker and its degree of clinical impact. As a result, an overall score of high, medium, or low is provided for each functional pillar. The score is represented by color-coded icons and informational graphics.

The specific biomarkers of concern that are utilized to establish the results for each functional pillar are listed in the Four Functional Pillars Biomarker Map.

Four Functional Pillars Biomarker Map

Infection	Inflammation	Insufficiency	Imbalance
any parasite present	Calprotectin	PE1	n-Butyrate
any pathogen present	EPX	Total Fecal Fats	Total SCFA
	Fecal IgA	Total Protein Products	Beta-glucuronidase
			<i>Lactobacillus</i>
			<i>Bifidobacterium</i>
			<i>E. coli</i>

Diversity and Abundance

It is now known that the human GI tract is home to more than 1000 species of microbial organisms, almost all of them bacteria. These organisms – collectively known as the microbiome – far outnumber the human cells in any individual and fulfill many metabolic functions.¹⁻³

It is becoming evident that many factors go into developing and maintaining what might be called a “healthy microbiome,” but this emerging area of biology is vastly complex. Indeed, prior to very modern analytic developments, there was no way to understand either the number or the functions of the tremendous population of organisms making up the human GI microbiota. Recent developments rely on DNA or RNA patterns, and by comparing detected sequences with libraries of known organisms, laboratories can now detect a tremendous number and variety of organisms.


One of the first fruits of this new technology is the discovery that, while the entire pool of possible members of the microbiota is large (> 1000 species), a much smaller number, 150 to 170 species, is found to predominate in any given individual.

Utilizing a molecular assay platform optimized for stool analysis, *GI Effects* assesses a set of clinically relevant 24 genera/species that map to 7 major phyla. *GI Effects* utilizes 16S rRNA gene polymerase chain reaction (PCR) amplification techniques. The improved method offers an enhanced DNA extraction method, updated primer sequences, and optimized thermal cycling for an enhanced reportable range. After detection of these organisms, a computer algorithm is used to map them into a graphic representations of diversity association and relative abundance.

Diversity Association is a proxy measure of gut biodiversity, which is defined as the number and abundance of distinct types of organisms present in the gut.⁴ The clinical utility of biodiversity of the gut microbiome is not yet fully defined, but in general a high diversity of GI organisms has been associated

with states of relatively good health, while low diversity has been associated with states of disease or chronic dysfunction in the scientific literature.

Relative Abundance represents the levels of selected phyla in an individual's microbiome and is represented relative to similar measures derived from a healthy cohort of individuals.



An algorithm is also used to graphically represent the microbial diversity and relative abundance of the commensal bacteria. The Diversity Association is a proxy measure of the diversity level of organisms in the gut. The impact of each commensal genera/species is weighted based on its association with high and low bacterial diversity. The Relative Abundance represents the size of each of the phyla as calculated by the assessed commensal organisms and is shown in comparison to the levels seen in a defined healthy cohort.

Measuring diversity association and relative abundance serve as global measures of an individual's gut health. Monitoring these markers over time provides insight into the impact of medications, supplements, diet, and lifestyle interventions on current and future health. Specific treatments, as directed by the results of the four functional pillars, would be expected to produce positive responses in the diversity and abundance of gut bacteria.

Biomarker Review

Modern technology allows the use of a growing number of biomarkers found in stool to supplement, and often supplant, more invasive and generally more expensive tests of GI function.

While the most obvious role of the human gastrointestinal (GI) tract is the incorporation of nutrients and energy from the diet, and the elimination of waste products and toxins, it is now clear that functions of the GI tract influence, not only GI health, but that of the entire human organism.

For these reasons, the results from the *GI Effects Comprehensive Profile* are reported using the “DIG” framework, which provides information on the three main categories of GI function:

- **Digestion and Absorption** reports on the effectiveness of GI breakdown and absorption of nutrients from ingested food. This category contains:
 - › Pancreatic Elastase-1 (PE1), a marker of exocrine pancreatic function
 - › Products of Protein Breakdown identifies bacterial fermentation of proteinaceous material
 - › Fecal Fat, a marker of fat maldigestion and malabsorption
- **Inflammation and Immunology** reports on the functioning of the inflammatory response and the secretory immune system in the GI tract. This category contains:
 - › Calprotectin, a marker of neutrophil activity and inflammation
 - › Eosinophil Protein X, a marker of eosinophil activity and inflammatory, allergic, and parasitic influences
 - › Fecal secretory IgA (fsIgA), a marker of secretory immune function, GI mucosal defense, and the maintenance of gut barrier function
- **Gastrointestinal Microbiome** reports on the status and function of the hundreds of microbial species (chiefly bacteria and fungi) that constitute the non-host living contents of the human GI tract. This category contains the following subcategories and their constituents:
 - › **Metabolic**, a series of biomarkers that indicate the microbial production of beneficial molecules including metabolites of incompletely-digested nutrients, as well as enzymes involved in the trafficking of biliary conjugate molecules. This subcategory contains:
 - Short-chain fatty acids (SCFAs), markers of colonic fermentation of carbohydrates into SCFAs required for colonocyte health and signaling
 - Beta-glucuronidase, an enzyme involved in carbohydrate digestion and lysis of bonds linking bioactive molecules to their glucuronide conjugates
 - › **Commensal Bacteria (Polymerase Chain Reaction – PCR)**, a listing of the levels of 24 major bacterial genera/species in 7 major phyla of commensal organisms, i.e., organisms known to exist symbiotically with humans, sharing and exchanging metabolic functions. Increasing evidence suggests that the makeup and metabolic functions of the commensal bacteria in the human microbiome are essential to maintenance of general homeostasis and health of the host organism.

- PCR can evaluate anaerobic targets and provides quantification of each target, giving a semi-quantitative result.
- A 25th biomarker is also provided in this section, the Firmicutes/Bacteroidetes Ratio (F/B Ratio). The F/B ratio provides an estimate of the predominance of two major phyla of commensal organisms, which has been associated with a number of metabolic disorders.
- › **Bacteriology (Culture)**, a listing of both commensal and additional bacteria grown in traditional culture media. This familiar suite of tests complements targeted commensal bacterial detection by PCR methods by identifying the presence of hundreds of additional bacterial species, including opportunistic/potentially pathogenic bacteria, which may be indicators of imbalance in the gut microbiome.
- › **Mycology (Culture)** detects fungal organisms using traditional culture techniques.
 - Cultured organisms, both bacteria and mycology, are definitively identified using MALDI-TOF (Matrix Assisted Laser Desorption Ionization Time-of-Flight) technology. The MALDI-TOF mass spectrometry platform utilized for the rapid identification of bacteria and yeast from pure cultures on the *GI Effects Comprehensive Profile* report relies on the most extensive FDA-cleared library of microbial targets available on the market.
 - Culture identifies all viable cultivable organisms, and is the accepted standard for assessment of aerobic organisms.
- › **Parasitology** detects intestinal parasites by means of two complementary techniques:
 - Ova and Parasites (O&P); microscopic examination using light microscopy
 - Parasitology Enzyme ImmunoAssay (EIA) tests, which can detect specific pathogenic organisms of interest; this test panel detects *Cryptosporidium* species, *Giardia lamblia*, and *Entamoeba histolytica*.
- › **Bacteria Sensitivity** lists cultured bacterial pathogens and their relative susceptibility to both prescription and natural antimicrobial agents.
- › **Mycology Sensitivity** lists cultured fungal pathogens and their relative susceptibility to both prescription and natural antimicrobial agents.

The *GI Effects Comprehensive Profile* report concludes with a section on Additional Results. These include time-honored characteristics of fecal specimens such as occult blood, color, and consistency, as well as the option to add on specific EIA testing for certain known pathogens, including *Helicobacter pylori* stool

antigen, *Campylobacter* species, *Clostridium difficile*, and Shiga toxin-producing strains of *E. coli*.

Organization of the Biomarker Review

In the following sections,

- Each **biomarker** is first identified and described
- The candidate **patient population** for the biomarker is described
- The **comparator** or existing gold standard for the biomarker is presented (if established), along with performance characteristics of the fecal biomarker when appropriate
- The **interpretation of the fecal biomarker** is discussed, including the significance of out-of-range results
- Desirable **outcomes and therapeutic recommendations** are discussed, indicating how the specific test might benefit patients in a variety of clinical settings




Graphical Representation and Color Coding for Biomarkers

In addition to a numeric result and a stated reference range for each result, all biomarker results on the *GI Effects* report are graphically represented and color coded in the context of a specific reference population by means of **quintile reporting** or **medical decision point reporting**.

Quintile reporting permits the clinician to recognize at a glance where each individual patient result falls compared to the distribution of results for the reference population, which is often (but not always) a representative sample of the entire population that the laboratory has tested for each biomarker. After rank ordering each individual test result from lowest to highest, the

reference population is divided into 5 equal groups, or **quintiles**, where each group represents 20% of the total count of individual results in the reference population.

The following example shows a patient's result (black diamond) that lies in the middle of the 4th quintile. The diamond placement below indicates that approximately 70% of all patients in the reference population for this biomarker had results that are lower than this patient's result.




The quintile reporting bar is denominated in percentile units and the color thresholds on the bar approximate standard deviation (SD) thresholds, assuming a normal distribution. In general, for a 2-tailed test as illustrated above, the green region includes plus or minus 1 SD from the population mean, or 68.2% of all results. The yellow areas include plus or minus 2 SD from the population mean, and encompass 95.4% of the distribution, and the red area represents the remainder of the population that falls outside of 2 SD in either direction.

By examining the quintile reporting bar for each biomarker in the context of the numerical result and reference range, the clinician can quickly identify results calling for closer clinical consideration. In conjunction with the patient's history and physical findings, biomarker results trending towards or outside of 2 SD may require additional evaluation.

Quintile reporting may also be useful when serial testing is used, to assess movement in either direction, and as a monitoring tool for effectiveness of interventions.

Some biomarkers have established threshold values associated with specific clinical conditions, histopathological findings, or recommended clinical interventions. These biomarkers are not reported using the quintile system but instead are reported on a colorimetric graphical bar denominated in the same laboratory measurement units as the biomarker result. The color-coded thresholds are defined by specific **Medical Decision Points**, or MDPs. For example, fecal calprotectin, PE1, and eosinophil protein X all have reference ranges based on a clinically characterized healthy reference population (i.e. not a symptomatic tested population) and cutoff points indicating a normal result, a borderline or weakly positive result, and an abnormal or strongly positive result.


An example MDP reporting bar for a normal calprotectin result is shown here:



This fecal calprotectin result of ≤ 50 mcg/g stool in a patient meeting Rome criteria for Irritable Bowel Syndrome (IBS) virtually excludes the likelihood of Inflammatory Bowel Disease; however, a fecal calprotectin result above 120 would be a strong indication for additional evaluation (perhaps including colonoscopy) to determine the specific etiology of the inflammatory process.

Since 88% of patients with GI symptoms have a fecal calprotectin result < 50 , a quintile graphical representation would restrict all clinically significant values (those above 50) to the right hand side of the reporting bar on the far side of the 5th quintile.

For the sake of illustration, consider how results reporting for a patient with an initial calprotectin value of 150 mcg/g – and subsequent value of 60 mcg/g – would appear in quintile and MDP representations.



Thus, the MDP colorimetric reporting bar facilitates visualization of a clinically significant improvement in fecal calprotectin levels (e.g. results from 150 mcg/g to 60 mcg/g) better than a quintile reporting bar (e.g. from 95th percentile to 90th percentile).⁵





Digestion and Absorption

For proper nutrition and gastrointestinal (GI) function, ingested nutrients must first be broken down (digested, in a biochemical sense), and the products of digestion must then be absorbed through a variety of physical and biochemical processes.

Biomarkers of digestion and absorption provide information about nutrient breakdown and entry into the circulation. They ultimately indicate how well the GI tract is performing its basic digestive functions.

The biomarkers are:

- **Pancreatic Elastase-1**, a marker of exocrine pancreatic function
- **Products of Protein Breakdown**, markers of undigested protein reaching the colon
- **Fecal Fat**, a marker of fat breakdown and absorption

In good health, digestion is accomplished in several steps. First by chewing and other physical processes, and then by the actions of stomach acid and a host of enzymes produced in the pancreas and small intestine, breaking down the three major components of food: complex carbohydrates (starches), proteins, and fats.

Absorption of the resulting products of digestion then occurs by several distinct processes. Damage to, or impairment of, any of the processes involved in digestion or absorption results in two main problems: inadequate net absorption of nutrients, producing absolute or relative nutrient deficiencies, and/or delivery of intact nutrients to the colon, where gut microbes may inappropriately digest or ferment nutrients. Such fermentation results in byproducts leading to excessive osmotic loads and gasses, leading to abdominal discomfort, diarrhea, flatulence, and other common symptoms.⁶

Maldigestion is defined as impaired breakdown of nutrients, and is often the result of inadequate or impaired digestive enzymes (or gastric acid production), while malabsorption refers to impairments in absorption of the normal end products of digestion.⁷

When faced with a patient experiencing sub-acute or chronic GI symptoms, it is the task of the clinician to discern which, if any, of these processes is occurring, and then, if possible, to identify one or more underlying, primary causes. Finally, in many cases, once a primary cause has been identified, a rational and usually simple course of therapy may be prescribed, with the goal of repairing the underlying pathological processes.

In many cases, fecal biomarker testing is useful in discerning whether maldigestion, malabsorption, or both, are present. Such testing is also helpful in identifying the underlying causes, for which treatment may be available.

Pancreatic Elastase-1 (PE1)

The Biomarker

- The exocrine portion of the pancreas (the cells and structures not related to endocrine functions, such as insulin production) secretes numerous digestive enzymes, among them PE1
- PE1 is a robust proteolytic enzyme that reaches the colon without itself being digested, is not greatly affected by increases or decreases in intestinal transit times, and is not affected by pancreatic enzyme replacement therapy
- The PE1 reference range was adopted from an FDA-approved kit

Biomarker Key Points

- Noninvasive biomarker of pancreatic exocrine (i.e., digestive) function
- Is not affected by supplemental pancreatic enzymes
- Reflects true pancreatic exocrine function⁸

Fecal PE1 testing can be used for initial determination of pancreatic exocrine insufficiency in suspected patients, as well as the monitoring of pancreatic exocrine function in patients under treatment.

Patient Populations of Interest

Patients in whom PE1 testing may be useful include those with:

- Unexplained diarrhea
- Weight loss
- Other symptoms of maldigestion
- Abdominal pain
 - Including symptoms meeting clinical criteria for irritable bowel syndrome (IBS)
- Low bone density

In addition, pancreatic exocrine insufficiency may occur secondary to:

- Chronic pancreatitis
- Diabetes
- Celiac disease
- Cystic fibrosis
- Inflammatory bowel disease (IBD)
- Excessive alcohol consumption
- Gallstones

Fecal PE1 testing can be used for **initial determination** of pancreatic exocrine insufficiency in suspected patients

There is also evidence that aging populations may exhibit a progressive loss of PE1, since pancreatic exocrine function may decrease with age.⁹⁻¹⁴

Comparator/Gold Standard Tests

PE1 has a strong correlation with the gold-standard test for pancreatic insufficiency, the secretin-caerulein test.^{8,15} (Caerulein is a synthetic analog for pancreozymin and stimulates pancreatic activity in a similar manner.)

Interpretation

Fecal PE1 Value (µg)	Interpretation ^{16,17}
> 200	Normal exocrine pancreatic function
100 to 200	Mild-to-moderate exocrine pancreatic insufficiency
< 100	Severe pancreatic insufficiency

Results are based on MDPs and are not represented by quintile values.

Notes on Interpretation

- Fecal PE1 testing may have reduced sensitivity for detecting mild pancreatic exocrine insufficiency in children.¹⁸
- Consumption of vegetarian or vegan diets, or other diets involving decreased meat intake, have been associated with reductions in fecal PE1.^{19,20}
- Pancreatic exocrine insufficiency occurs in about 50% of type 1 diabetics, and in about 33% of type 2 diabetics.²¹
- Chronic pancreatitis patients may have compromised antioxidant systems.⁶

Outcomes and General Therapeutic Considerations

Patients with PE1 results suggestive of exocrine pancreatic insufficiency should undergo further investigation to determine the underlying causes of their dysfunction, as shown in the following table:

Additional Testing	Rationale
Evaluation of fecal fats ²²	Excess fecal fat may be due to: <ul style="list-style-type: none"> • Lack of bile acids (due to liver damage, hypolipidemic drugs, or impaired gallbladder function) • Celiac disease • Small bowel bacterial overgrowth • Other conditions and medications (e.g., Orlistat)
Full nutritional assessment	Defective exocrine pancreatic function may be associated with: <ul style="list-style-type: none"> • Abnormal blood lipids • Low levels of minerals (magnesium, zinc, selenium, and calcium) • Low levels of fat soluble vitamins (A, D, E, and K)²³

Supporting the Patient with Evidence of Pancreatic Exocrine Insufficiency

Certain lifestyle, medication, and supplement interventions may be appropriate for patients with abnormal fecal PE1 results suggestive of pancreatic exocrine insufficiency.

Lifestyle Support²³

- Small, frequent meals (better absorbed)
- Reduce alcohol consumption
- Smoking cessation

Medication/Supplement Support^{22,24-26}

- Support patients with pancreatic exocrine insufficiency by pancreatic enzyme replacement therapy (PERT) at doses appropriate for degree of insufficiency and based on symptom improvement; in some conditions PE1 levels normalize as underlying disorders improve²⁷ (improved PE1 levels reflect functional improvements, not supplementary enzymes)





Products of Protein Breakdown

The Biomarker

When proteins or their digestion products (oligopeptides and amino acids) reach the distal colon, they are fermented by colonic organisms (proteolytic fermentation) into a group of compounds including the characteristic short-chain fatty acids (SCFAs), isovalerate, isobutyrate, and valerate.

Biomarker Key Points

- Normal protein digestion and absorption is relatively complete in stomach and small intestine
- Healthy colonic contents therefore include only small amounts of protein-derived SCFAs
- Protein fermentation can yield a diversity of end products, including SCFAs, amines, phenols, indoles, thiols, sulfur compounds and branched-chain fatty acids²⁸
 - Though many have shown toxic properties in vitro and in animal models, the relationship between gut health and protein fermentation in humans has not been thoroughly investigated²⁹⁻³¹
- Primary colonic SCFAs from protein breakdown are valerate, isovalerate, and isobutyrate^{16,17}
- The result on the *GI Effects* report reflects a combined total of valerate, isovalerate, and isobutyrate measurements

Patient Populations of Interest

Patients with protein maldigestion, or those with abnormally large amounts of protein presented to the distal colon, may demonstrate increased products of colonic protein breakdown in the stool.

Increased fecal presence of products of protein breakdown may be present in patients with:

- Hypochlorhydria (diminished hydrochloric acid secretion in the stomach), which is associated with:³²
 - › Advanced age: in roughly 30% of elderly patients, gastric acid secretion is diminished
 - › Use of acid-blocking medications or dietary supplements that produce too high a gastric pH to allow for complete protein digestion in the stomach³³
 - › Food reactions; elevated gastric pH (less acidic) has been associated with increased risk for food reactions, possibly from hindering protein breakdown³⁴⁻³⁶

- Pancreatic exocrine insufficiency or pancreatitis; insufficient pancreatic proteases leave improperly-digested protein fragments that reach the colon^{37,38}
- Excessive protein intake³⁹
- Gastrointestinal bleeding or irritation, mucosal desquamation, and bacterial overgrowth; these conditions result in excessive self-derived proteins in the intestinal lumen⁴⁰

Comparator/Gold Standard Tests

Currently there is no gold standard assessment for fecal products of protein breakdown. Products of Protein Breakdown are utilized as a contributory diagnostic tool. Testing for fecal nitrogen may indicate the presence of protein malabsorption, but fecal nitrogen is difficult to measure and is not in widespread clinical use.⁴¹

Interpretation

The result for products of protein breakdown reflects the sum of fecal valerate, isovalerate, and isobutyrate. On the *GI Effects* test report, the value is shown against a background representing population quintiles, as described in the introduction.

Outcomes and General Therapeutic Considerations

Patients with elevated products of protein breakdown should be evaluated for common causes of insufficient protein digestion and/or excessive protein presenting to the colon.

Source of Elevated Colonic Products of Protein Breakdown	Possible Causes	Therapeutic Response
Insufficient Protein Digestion ^{33,34}	Hypochlorhydria	<ul style="list-style-type: none"> • Reduce acid-blocking medications • Add betaine HCl
	Pancreatic exocrine insufficiency	<ul style="list-style-type: none"> • Evaluate fecal PE1
Excessive delivery of protein to colon ⁴²	High-protein diet	Review protein/carbohydrate intake
	GI irritation/inflammation, bleeding, bacterial overgrowth	Additional testing, e.g., fecal calprotectin, fecal eosinophil protein X (EPX), fecal occult blood, stool culture for beneficial bacteria

Total Fecal Fats

The Biomarker

Under normal conditions, the bulk of dietary fat is digested and absorbed in the small intestine, leaving only small amounts for delivery to the colon and fecal stream. Fecal fat measurements determine the amount of fat in stool, and may therefore identify fat maldigestion, malabsorption, or steatorrhea.

Biomarker Key Points

- The test is a fecal fat extraction method that results in a quantitative value.
- Fecal fat extraction methods have been found to correlate with degree of fat malabsorption.⁴³
- Total fecal fat is made up of long-chain fatty acids (LCFAs), cholesterol, triglycerides, and phospholipids.

Patient Populations of Interest

Fecal fats should be measured in any patient for whom steatorrhea (passage of pale, bulky, and malodorous stools) may be a symptom of underlying digestive or non-digestive disorders.

Symptoms suggesting evaluation of fecal fat as a means of detecting root causes include:

- Fatigue
- Unexplained anemia
- Nutrient deficiencies
- Unintended weight loss

Comparator/Gold Standard Tests

The 3-day stool collection with total fecal fat determination is the gold standard test for fecal fat. This test is unwieldy and unpleasant for patients and lab personnel. The total fecal fat extraction on a single specimen provides a quantitative value to identify patients that may benefit from the more in-depth 3-day test. Limited research has found extraction methods to correlate with the gold-standard.⁴³

Interpretation

Total fecal fat is the sum of fecal triglycerides, long-chain fatty acids, cholesterol, and phospholipids. Fecal fats are reported using the quintile system.

Outcomes and General Therapeutic Considerations

Fecal fat may be elevated in situations of *fat maldigestion*, such as:

- Pancreatic exocrine insufficiency (inadequate lipase production or delivery)
 - › Causes include chronic pancreatitis and cystic fibrosis
- Bile salt insufficiency (inadequate solubilization of fats for digestion)
 - › Causes include liver damage, hypolipidemic drugs, impaired gallbladder function
- Hypochlorhydria (inadequate stomach acid)
 - › Causes include aging and gastric acid-lowering drugs
- Small intestinal bacterial overgrowth and resulting acidic small-intestinal pH (impairment of small intestinal digestive enzymes)^{15,30}
- Use of medications designed to impair intestinal lipase activity (Orlistat, Xenical, Alli), or use of synthetic fat-like products indigestible by normal lipase (Olestra)⁴⁴⁻⁴⁷
- Elevated fecal fat may be associated with deficiencies in fat-soluble nutrients, so consider nutritional assessment of essential fatty acids, fat-soluble vitamins, and minerals.^{48,49}

Fecal fat may also be elevated in situations of *fat malabsorption*, such as:

- Intestinal dysbiosis
- Intestinal parasites
- Gastric bypass, ileal resection, or other surgeries that limit absorptive surface area

Finally, fecal fat may be elevated in patients with:

- Irritable bowel syndrome (often as a symptom of pancreatic exocrine insufficiency)
- Inflammatory bowel disease
- Food intolerances⁵⁰
- Celiac disease
- Excessive alcohol intake
- Chronic use of non-steroidal anti-inflammatory drugs (NSAID)⁵¹



Supporting the Patient with Elevated Fecal Fat Levels

Depending on root causes, patients with elevated fecal fat levels can be supported as follows: ^{6,45-47,52-54}

Suspected Cause of Elevated Fecal Fat	Support Measures	Rationale
Pancreatic exocrine insufficiency	Supplementary plant or pancreatic digestive enzymes	Lipases increase fat digestion
Disorders of bile formation/transport (cholestasis)	<ul style="list-style-type: none"> Bile salts or cholagogues, taurine or glycine Diet changes 	Enhance intestinal fat solubilization
Lipase inhibitors (orlistat, Xenical, Alli) or synthetic fat consumption (Olestra)	Discontinue these products	Permit normal fat digestion

Inflammation and Immunology

Interactions between the immune system and the GI tract are being recognized as of growing importance, not only in GI physiology and pathophysiology, but also in their influences on systemic health and disease.

Biomarkers of GI inflammation and immunology provide information about the GI tract's interactions with, and responses to, the outside world. They indicate how well the GI tract is maintaining its role as a barrier, as well as whether the GI tract is undergoing pathological responses to external or internal challenges.

The biomarkers are:

- **Calprotectin**, a marker of neutrophil-driven inflammation
- **Eosinophil Protein X**, a marker of eosinophil-driven inflammation and allergic response
- **Fecal Secretory IgA**, a marker of gut secretory immunity and barrier function

Calprotectin The Biomarker

Calprotectin is a protein produced in abundance by neutrophils, the ubiquitous immune system “first responders.” When neutrophils accumulate at sites of inflammation, they release increased amounts of calprotectin in a way that closely correlates with findings on endoscopy and histology, and are thus useful in quantifying the degree of intestinal inflammation. ^{55,56}

This property makes calprotectin useful for differentiating inflammatory from non-inflammatory disease processes, e.g., distinguishing irritable bowel syndrome (IBS) from inflammatory bowel disease (IBD). ⁵⁷

Biomarker Key Points

- Calprotectin is described in the literature as a useful non-invasive screening tool for identifying which patients may benefit from endoscopy for suspected IBD
- Calprotectin is used in diagnosing IBD (Crohn's disease and ulcerative colitis) and for quantifying degree of inflammation
 - › Thus calprotectin may be useful for monitoring treatment and assessing for relapse in patients with known IBD
- Calprotectin is FDA-cleared to differentiate IBS from IBD ⁵⁷

Patient Populations of Interest

Patients with symptoms consistent with IBS should have fecal calprotectin testing done as a means of ruling out significant inflammation; those with positive Rome criteria and normal calprotectin (≤ 50 mcg/g) have virtually no chance of having IBD. ^{58,59}

Rome III Diagnostic Criteria for IBS ^{23,60*}

Recurrent abdominal pain or discomfort** at least 3 days/month in the last 3 months associated with two or more of the following:

1. Improvement with defecation
2. Onset associated with a change in frequency of stool
3. Onset associated with a change in form (appearance) of stool

* Criteria fulfilled for the last 3 months with symptom onset at least 6 months prior to diagnosis

** “Discomfort” means an uncomfortable sensation not described as pain

Studies suggest that a person with positive Rome criteria and a normal calprotectin (≤ 50 μ g/g) has virtually no chance of having IBD

Comparator/Gold Standard Tests

The gold-standard comparator test for determining presence and degree of intestinal inflammation is endoscopy with biopsy and histology; fecal calprotectin correlates closely with this approach.

Interpretation

The expected values for fecal calprotectin are shown here:

Calprotectin Concentration (micrograms/g stool)	Interpretation	Follow-up
≤ 50	Normal (no active GI Inflammation)	None
> 50 to 120	Borderline, suggestive of low-grade inflammation	Re-evaluate in 4–6 weeks
> 120	Abnormal	Determine source of inflammation and repeat as clinically indicated
> 250	Associated with high risk of clinical relapse in IBD ⁶¹	Adjust therapy accordingly

On the *GI Effects* report, calprotectin levels are represented as normal, borderline, or high, and are not represented by quintile values.

- False negatives for calprotectin may be seen in patients with severe immune compromise, who are not able to mobilize neutrophils sufficient to raise levels of calprotectin in the intestine
- Calprotectin elevations are seen in conditions other than IBD, e.g., malignancy, infection; therefore, a primary diagnosis of IBD cannot be established solely on the basis of a positive fecal calprotectin result
- Elevated calprotectin may also result from chronic NSAID use; evaluate such use in all patients with elevated calprotectin

Outcomes and General Therapeutic Considerations

Calprotectin is a simple, reliable, and non-invasive test that is useful in:

- Selecting patients with abdominal symptoms who may require further diagnostic procedures⁶²
- Aiding in distinguishing between IBD and IBS⁶³
- Selecting/screening patients for endoscopy, especially children in whom general anesthesia might be required for invasive study⁵⁵
- Determining disease activity and risk of relapse in IBD⁶⁴
 - › Therapy can then be initiated before inflammation reaches critical intensity
- Monitoring IBD treatment response and determining when a full clinical remission has been achieved⁶⁵
- Evaluating efficacy in trials of new treatments for IBD

Supporting the Patient with Elevated Fecal Calprotectin Levels

- Because inflammatory, infectious, or neoplastic processes may result in an elevated calprotectin level, the cause of a value > 120 mcg/g warrants further investigation – including endoscopy or radiography – based on clinical correlation
 - › Assessments to uncover causes of bowel inflammation are found at the end of this section on Inflammation and Immunology

Eosinophil Protein X (EPX) The Biomarker

Many inflammatory and neoplastic processes in the gut involve increased activity of eosinophils, white blood cells that normally reside in the lamina propria (connective tissue layer) of the intestinal wall.

When the lamina propria is damaged, eosinophils migrate into the gut lumen, where they degranulate to release a variety of proteins with cytotoxic properties, which contribute to ongoing inflammation and tissue destruction.⁶⁶ One such protein is eosinophil protein X (EPX), which can be measured in fecal matter.⁶⁷

Biomarker Key Points

- Fecal EPX offers the practitioner a noninvasive alternative to the invasive gold standard, allowing for better differential diagnosis
- EPX is considered the superior cationic protein for assessment of eosinophil function, because it most accurately reflects the degree of mucosal damage⁶⁸
- Baseline EPX levels offer a way to determine and monitor GI inflammation associated with food allergy
 - › Significant reduction in EPX after 3 months on an elimination diet has been demonstrated⁶⁹
 - › Return of EPX to normal levels (< 2.0 mcg/g) can be used to indicate clinical efficacy of elimination diets⁶⁹ or clinical remission of IBD⁶⁸
- Systemic corticosteroids can reduce circulating levels of EPX⁷⁰

Patient Populations of Interest

Clinically, elevations in EPX indicate the presence of an IgE-mediated inflammatory process. Patients at risk for having such processes, and in whom fecal EPX may be helpful in making diagnosis and treatment plans, include those with:

- Complaints of GI symptoms related to food intake (possible food allergy)
 - › Also for monitoring results of elimination diets or other interventions in patients with known food allergy
- Concerns about possible parasitic/worm infections (recent travelers with GI symptoms)
- IBD in need of non-invasive monitoring for disease activity and treatment monitoring (consider calprotectin as well)

Because food allergy may present with symptoms similar to those seen in IBS, also consider fecal EPX testing in patients with symptoms consistent with IBS.

Comparator/Gold Standard Tests

The gold standard test for quantifying eosinophils in the gut is whole gut lavage. This is an invasive procedure requiring endoscopy, and has limited utility in the office setting. Clinical research indicates significant correlation between whole gut lavage fluids and eosinophil mediators in stool, such as fecal EPX.⁷¹

Interpretation

The normal reference range for EPX is < 7.0 mcg/g of stool. Levels above 7.0 suggest increased eosinophil activity in the gut lumen, suggestive of allergic or inflammatory reactions.

- Baseline EPX levels may be used to determine inflammation associated with food allergy
 - › Subsequent testing can be used to monitor results of dietary changes
- Elevations of EPX correlate with disease activity in ulcerative colitis and Crohn's disease⁶⁸
- Serial testing of EPX offers a non-invasive means of evaluating disease activity and for predicting relapses in patients with IBD⁷²

Outcomes and General Therapeutic Considerations

EPX levels above 2.0 mcg/g of stool are suggestive of elevated eosinophil activity, which may be associated with the following:^{68,71,73}

- IgE-mediated food allergy
- Intestinal parasitic infection
- IBD

Less commonly, fecal EPX may be elevated in:^{66,71,72,74,75}

- Atopic dermatitis
- Gastroesophageal reflux disease (GERD)
- Collagenous colitis
- Allergic colitis
- Excessive alcohol intake
- Chronic diarrhea
- Protein-sensitive enteropathy
- Gastrointestinal cancer
- Eosinophilic gastroenteritis (rare)

When evaluating a differential diagnosis that includes bowel cancer, calprotectin may be useful as an adjunctive marker along with radiologic and/or endoscopic evaluation

Conversely, in patients with known food allergy or IBD, acquisition of a normal EPX (< 2.0 mcg/g stool) can indicate the efficacy of treatment, such as elimination diet or remission of IBD.

Supporting the Patient with Elevated Fecal EPX

To further delineate the root cause(s) of an elevated fecal EPX result, consider additional testing to uncover causes of bowel inflammation, which are found at the end of this section on Inflammation and Immunology.



Fecal Secretory IgA (sIgA) The Biomarker

Secretory IgA, or sIgA, is a class of antibodies produced by and secreted from mucosal surfaces, especially the GI and respiratory tracts. In the gastrointestinal epithelium, sIgA is the first line of defense against the entry of enteric toxins and pathogenic organisms from the colon. Colonic sIgA is closely involved in maintenance of the gut epithelial barrier, and in the development of immune tolerance of normal, beneficial commensal gut organisms, as well as of common molecular epitopes found in foods.^{76,77}

Measurement of sIgA in fecal material, or fecal sIgA, therefore, may provide longitudinal information about the status of the intestinal epithelial barrier and its function in both immune exclusion (prevention of pathological material and organisms from gaining entry into the general circulation) and in immune inclusion (delivery of commensal bacteria and their products to the gut and systemic immune system for recognition and the development of tolerance, which spares these beneficial organisms from destruction by the immune system).^{78,79} To date, there is little support from scientific data for the use of sIgA as a stand-alone diagnostic study; rather, measurement of sIgA has been most valuable in determining the outcomes of selected interventions that modulate GI barrier functions, including both disturbances of immune exclusion (in which excessive foreign matter is allowed to enter the organism, resulting in colonization by pathogens) and immune inclusion (in which immune tolerance fails or is weakened, resulting in allergic manifestations).⁷⁹

Biomarker Key Points

- Fecal sIgA in the gut:
 - › Regulates the balance of commensal (beneficial) bacteria^{80,81}
 - › Prevents colonization by pathogens⁸²
 - › Promotes tolerance of commensal organisms^{76,77}
 - › Maintains GI barrier function⁷⁸
 - › Promotes formation of normal biofilm containing beneficial organisms in the gut⁸³

Patient Populations of Interest

Patients with known or suspected disruptions of the GI epithelial barrier may suffer from manifestations of “leaky gut” such as excessive bacterial translocation (failure of immune exclusion) or food allergies (failures of immune inclusion and tolerance development).⁸³⁻⁹⁰

Management of such patients may involve treatment with probiotics for restoration of a healthier colonic flora, or with exclusion diets in attempts to identify and eliminate sources of

food allergies/intolerances. Measurement of sIgA is a common means of monitoring response to therapy in studies of such interventions, and clinicians may find it useful to determine sIgA levels prior to and during the course of treatment.⁹¹⁻⁹⁴

The sIgA study is not recommended for use as a primary diagnostic test in patients suspected of having congenital immunoglobulin deficiencies such as selective IgA deficiency; such immune deficiencies should be evaluated by serological tests based on symptoms and family history.⁹⁵

Comparator/Gold Standard Tests

There are no standard comparator tests for secretory IgA production in the gut.

Interpretation

Fecal sIgA results are reported on the *GI Effects Comprehensive Profile* using the quintile reporting system.

- There are no validated correlations between sIgA levels and specific disease states
- Fecal sIgA is most commonly used as an outcomes measure in clinical trials
- Determination of a pre-intervention sIgA level, followed by serial post-intervention levels, may clarify responses to food elimination diets, pre- and probiotic supplementation regimens, and similar therapeutic maneuvers
- Increased sIgA levels may indicate a normal, transient, immunological response to intestinal viral or bacterial pathogens⁹⁶
 - › Monitoring the level may aid in identifying resolution of issue

Outcomes and General Therapeutic Considerations

Observation of the sIgA levels following therapeutic interventions may help to determine whether and how much a patient is responding to such therapies. Fecal sIgA levels should not be used in isolation to make diagnostic or therapeutic decisions, but rather should form part of the picture of gut barrier function in the context of additional clinical data.



Lactoferrin (Add-on) The Biomarker

Lactoferrin is an iron-binding glycoprotein secreted by most mucosal membranes; it is a major granular component of neutrophils (white blood cells). Liberated from the neutrophils in response to inflammation, lactoferrin binds to iron, impeding microbial growth and facilitating generation of hydroxyl radicals.⁹⁷

Biomarker Key Points⁹⁸

- Expressed by surface epithelial cells and found in most exocrine secretions, including breast milk, tears, nasal secretions, saliva, intestinal mucus, and genital secretions
- Lactoferrin is a multifunctional protein with antibacterial and immune modulatory activities and is a component of the first line of host defense
- Its expression is upregulated in response to inflammatory stimuli
- In the gastrointestinal tract, lactoferrin serves as a non-specific marker of inflammation

Patient Populations of Interest

Potential indications for testing are patients exhibiting gastrointestinal symptoms with suspected inflammation. Although clinical cut-offs have as yet to be determined definitively, lactoferrin may be useful in assisting the clinician in (1) identifying which patients may need further evaluation for inflammatory bowel disease⁹⁹⁻¹⁰¹ and (2) differentiating between IBD and non-inflammatory irritable bowel syndrome (IBS).¹⁰²

Comparator/Gold Standard Tests

The gold-standard comparator test for determining presence and degree of intestinal inflammation is endoscopy with biopsy and histology. Fecal lactoferrin has been correlated to histological findings.¹⁰³

Interpretation

A positive lactoferrin test generally indicates inflammation of the intestinal mucosa.

Outcomes and General Therapeutic Considerations

Lactoferrin is a non-invasive screening test and can assist the physician in stratifying symptomatic patients who require further evaluation; subsequent calprotectin testing can provide additional useful diagnostic information and assist in triage for endoscopic referral.

Supporting the Patient with Abnormal Lactoferrin Results

To further delineate the root cause(s) of an abnormal lactoferrin result, consider additional testing to uncover causes of bowel inflammation, which are found immediately below.

Supporting the Patient with Gastrointestinal Inflammation and Immune Reactions

Regardless of the specific findings, the support for patients with evidence of GI inflammation and/or immune reactions involves the following steps:

- Eradicate known pathogens or other infectious agents¹⁰⁴
- Consider supporting commensal bacteria with probiotic supplements and dietary changes^{92,105-107}
- Consider intestinal mucosal and anti-inflammatory support: appropriate nutrients, dietary changes, and botanicals¹⁰⁸⁻¹¹⁵
- Rule out food sensitivities or allergies; consider elimination diet and/or IgG and IgE food sensitivity testing¹¹⁶⁻¹¹⁸
- Support immune status with supplemental whey protein or increased fiber

Consider further evaluation of underlying causes as shown in the following table:

Tests for Discerning Underlying Causes of Bowel Inflammation

Intestinal Permeability (IP) Assessment	IP is a noninvasive assessment of impaired permeability or leaky gut, which can affect barrier and immune function. Chronic irritation in the gut lining can lead to maldigestion and malabsorption.
Food Antibody Assessment	Assessment of immune-mediated responses to specific foods, both IgE and IgG classes.
Celiac Panel	Celiac disease is an autoimmune response to the gluten protein, gliadin, and to proteins with similar structure; damage occurs to the villi in the small intestine. The panel includes highly sensitive markers for identifying celiac disease.
ImmunoGenomic™ Profile	Evaluates genetic variations that modulate immune and inflammatory activity; these polymorphisms can result in stimulation of mechanisms that lead to chronic, overactive inflammatory responses.

Gastrointestinal Microbiome

The gastrointestinal microbiome is emerging as an exciting and powerful area for managing not only GI health, but that of the entire organism. It is estimated that the GI microbiome contains at least 10 times as many individual cells as are found in the somatic cells of human beings, with at least 1000 times the amount of genetic information.

Biomarkers of the GI Microbiome provide information about the health, function, and diversity of the trillions of microbial cells in the GI tract. They indicate how well the microbiome is performing its shared metabolic functions with the human host.

The biomarkers are:

- **Metabolic indicators**, which demonstrate specific and vital metabolic functions performed by the microbiota
- **Commensal Bacteria**, which demonstrate the composition, diversity, and relative abundance of gut organisms, all of which are being linked to general health
- **Bacterial and mycological culture**, which demonstrate presence of specific beneficial and pathological organisms
- **Parasitology**, which demonstrates presence of parasites

Clinically significant imbalances in the GI microbiome (dysbiosis) have been associated with:¹¹⁹⁻¹²¹

- Irritable bowel syndrome (IBS)
- Inflammation
- Inflammatory bowel disease (IBD)
- Immune Modulation
- Metabolic Disorders
- Body weight & fat distribution
- Insulin sensitivity/type 2 diabetes
- Autism
- Other acute and chronic disorders^{23,61,122-126}

There are multiple ways of assessing the GI microbiome and its impact on human health.

Humans have a conserved set of gut microbiota that is generally shared by most individuals, though each person has their own distinct and highly variable microbiota. The commensal gut microbiota interacts extensively with the host, influencing multiple metabolic and physiological functions.

Metabolomic methods involve measuring and interpreting the metabolic products of bacterial activity in the colon; these include determinations of fecal short-chain fatty acids and of enzymes elaborated by gut organisms, such as beta-glucuronidase.

Taxonomic methods are rapidly maturing, and include traditional bacterial and fungal culture and microscopy, as well as newer methods based on the microbial genomes of individual bacterial groups, such as polymerase chain reaction (PCR).

Metabolic Products: Short-Chain Fatty Acids The Biomarker

Commensal gut bacteria anaerobically ferment resistant starch and dietary fiber (including prebiotics) to produce the beneficial short-chain fatty acids (SCFAs) acetate, propionate, and butyrate.¹²⁷ In particular, n-butyrate is the obligate fuel source for colonocytes, and inadequate levels are associated with disordered colonic health.^{128,129}

Studies show distinct differences between healthy populations and those with GI disorders in terms of gut microbiota composition and SCFA production and distribution.

In the gut, short-chain fatty acids:

- Maintain intestinal barrier function
- Provide fuel for colonocytes (n-butyrate)
- Regulate colonic absorption of water and electrolytes
- Salvage unabsorbed carbohydrates
- Support commensal bacteria

SCFA levels are influenced by many factors, including:

- Diet composition
- Fecal ammonia content (indicative of excessive undigested protein)
- Obesity
- Environment

Biomarker Key Points

- Fecal SCFA concentrations are a metabolomic indicator of the health of the GI microbiome
 - › Low concentrations suggest either dysbiosis (abnormal levels or function of gut bacteria) or inadequate substrate for commensal organisms to ferment
- While optimal levels have not yet been identified, higher fecal concentrations of SCFAs have been associated with decreased GI disease in epidemiologic studies



Patient Populations of Interest

Studies have demonstrated distinct differences between healthy populations and those with gastrointestinal disorders when comparing the composition of gut microbiota and SCFA distribution. Lower levels of SCFAs are associated with colorectal diseases in population studies. Altered SCFA production has been seen in IBD and IBS. Fecal SCFA testing is a non-invasive contributory diagnostic tool in evaluating such patients.

Comparator/Gold Standard Tests

There is currently no gold-standard test for assessing SCFA production. SCFAs can also be measured in the blood. All three major SCFAs [acetate, butyrate, and propionate] are present in portal blood at concentrations several times greater than peripheral venous blood indicating the gut as a major source of these fatty acids.¹³⁰ Because SCFAs are rapidly utilized in the body, it is problematic to correlate blood levels with fecal SCFAs.

Interpretation

Optimal levels have not yet been determined for fecal SCFAs, however, in general, higher levels are considered beneficial.

The relative concentrations of n-butyrate, acetate, and propionate are also reported as percentages of total SCFAs. All results are reported and graphically displayed using the quintile reporting system.

- Low SCFA production is associated with:
 - › Decreased carbohydrate or fiber intake¹²⁷
 - › Low levels of fecal anaerobic or commensal bacteria (dysbiosis)^{62,131,132}
 - › Dysbiosis has been associated in research with inflammatory processes (e.g., IBD),¹³³ and functional bowel disorders (e.g., IBS)¹³⁴⁻¹³⁶

Outcomes and General Therapeutic Considerations

Low fecal SCFA levels typically indicate disordered metabolic processes in the colonic commensal community, e.g., inadequate amounts of beneficial bacteria or inadequate substrate for those bacteria to produce their beneficial metabolic products, SCFAs.

Diets high in fiber and resistant starch, and relatively low in protein content, may increase SCFA production with resulting reduced risk of colorectal diseases and lowering of cholesterol and blood sugar.¹³⁷⁻¹⁴⁵

Supporting the Patient with Low Fecal SCFA Concentrations

Support of patients with low fecal SCFAs is centered on increasing metabolic substrates for beneficial SCFA-producing organisms:

- Prebiotic supplementation
- Increased dietary carbohydrate and fiber intake
- Increased consumption of resistant starch, which is known to increase levels of fecal butyrate¹⁴⁶
- Evaluate and treat abnormalities of commensal gut bacteria

Beta-glucuronidase The Biomarker

Beta-glucuronidase is an enzyme that breaks down complex carbohydrates. Additionally, it acts to deconjugate glucuronide molecules from a variety of toxins, carcinogens, hormones, and drugs, which are naturally glucuronidated in the liver to facilitate biliary excretion. Deconjugation of these molecules in the gut permits their reabsorption via enterohepatic recirculation, producing higher than desired blood levels of potentially harmful compounds.

Beta-glucuronidase activity must be sufficient to permit deconjugation and absorption of desirable molecules, while remaining low enough to prevent widespread deconjugation and subsequent reabsorption of toxins and other undesirable molecules.

Additionally, many beneficial nutrients are ingested as the glucuronide conjugate of the active molecule, which must be deconjugated in order for the beneficial molecule (the “aglycone”) to be absorbed. Such nutrients include lignans, flavonoids, ceramides, and glycyrrhetic acid.

Thus, a proper balance of glucuronidase in the gut lumen is essential. Beta-glucuronidase is inducible in colonocytes, but it is also produced by anaerobic gut bacteria (particularly *E. coli*, but also *Peptostreptococcus*, *Bacteroides*, and *Clostridium*).

Biomarker Key Points

- Limited research suggests an association between elevated fecal beta-glucuronidase and colon cancer risk^{64,147,148}
- Low fecal beta-glucuronidase may also represent a problem, because the enzyme is needed to release the active aglycone forms of many dietary phytonutrients

Patient Populations of Interest

Evaluating beta-glucuronidase may be of interest to clinicians interested in evaluating substances that require deconjugation of glucuronide molecules, such as hormones, vitamin D, toxins, and phytonutrients.

Comparator/Gold Standard Tests

There is currently no gold-standard test.

Interpretation

Results are reported using the quintile reporting system.

Abnormally high levels of this biomarker warrant further investigation; abnormally low levels may diminish the bioavailability of many phytonutrients.

Outcomes and General Therapeutic Considerations

Further evaluation of patients with elevated fecal beta-glucuronidase includes consideration of exposure to and intake of toxins, hormones, and drugs.

For patients with persistently low fecal beta-glucuronidase, consider genomic testing and detoxification profiles to assess genetic markers and functioning of the glucuronidation pathway.

Supporting the Patient with Elevated fecal beta-glucuronidase:

For patients with **elevated** fecal beta-glucuronidase:


- The following supplements may be helpful:
 - › Calcium-D-glucarate
 - › Milk thistle
 - › Probiotics (*Lactobacilli* and *Bifidobacteria*)
- The following dietary management may be helpful:
 - › Increased consumption of vegetables and insoluble fiber

Commensal Bacteria

The Biomarkers

Commensal bacteria populate the human GI tract, especially the colon.

More than 95% of commensal gut organisms are anaerobic. Because these organisms are therefore difficult to recover by traditional (aerobic) culture techniques, there has been a major shift in research towards molecular DNA techniques, which are now considered the standard for anaerobic bacteria assessment in research. These techniques permit identification and quantification of multiple organisms with a single specimen.



The Polymerase Chain Reaction (PCR) methodology is capable of identifying and quantifying organisms. PCR probes can identify bacterial populations at any level of taxonomy, as broadly as phylum and as narrowly as species. This ability permits analysis of the gut microbiome at any desired degree of complexity.

Such analyses are of growing importance in human health and disease. The commensal gut microbiome has a diverse set of direct interactions with the human host. These interactions influence multiple metabolic and physiological functions, including:

- Producing short-chain fatty acids that nourish colonocytes and modulate gut physiology
- Modulating the systemic and intestinal immune systems
- Modulating GI hormone production
- Maintaining gut barrier function and motility
- Modulating oxidative responses





- Producing vitamins (e.g., biotin, vitamin K)
- Metabolizing xenobiotics and phytochemicals
- Preventing colonization by potential pathogens

In addition, the normal gut microbiome as a whole is important in preventing colonization by pathogens, a concept known as “colonization resistance.”^{149,150}

Clinically significant disruption or impairment of the healthy microbiome is referred to as “dysbiosis,” although there is as yet no standard definition of this term, nor is there an identified “ideal” or “normal” microbial community.¹⁵¹⁻¹⁶³ Some researchers suggest that dysbiosis may occur in three primary types: transient, disrupting, or persistent.¹³³

Biomarker Key Points

Assessment of an individual’s commensal bacterial population may provide clinical utility by permitting:

- Evaluation of 24 key bacterial groups/species that are associated with GI and systemic health
- Detection and characterization of states of imbalance between these groups
- Monitoring of changes following clinical interventions (e.g., antibiotic, probiotic, or dietary therapies)
- Tailoring of treatments to individual needs
- Identification of possible origins of metabolic disturbances that lie in alterations of the gut microbiome
- Detection of dysbiosis
 - › Commensal bacteria should be diverse and balanced in terms of population levels
 - › A low or imbalanced population of commensal bacteria may be an indicator of dysbiosis¹³³

Patient Populations of Interest

Because of the large numbers of conditions that are associated with imbalances in the composition of the gut microbiome, many patient populations may benefit from evaluation of GI commensal organisms. These include patients suffering from:

- Irritable bowel syndrome (IBS)
- Inflammation
- Inflammatory bowel disease (IBD)
- Immune Modulation
- Metabolic Disorders

- Body weight & fat distribution
- Insulin sensitivity/type 2 diabetes
- Autism
- Other acute and chronic disorders^{23,61,122-126}
- Conditions associated with “leaky gut” or impaired intestinal barrier function

Comparator/Gold Standard Tests

No gold standard exists for evaluation of the commensal bacterial community. DNA-based techniques such as PCR are now considered the standard for anaerobic bacteria assessment in research and such techniques are now being applied in the clinical realm. These techniques allow for vastly increased detection and quantitation of large numbers of organisms across the entire phylogenetic spectrum.

Interpretation


The commensal bacteria PCR testing panel produces results for 24 important bacterial groups/species known to be associated with both GI and systemic health.

- Organisms are reported according to phylum, and sub-categorized by specific genus and, in most cases, by species
- The abbreviation “spp.” is used to indicate multiple species within a genus
- Commensal bacteria PCR results are reported as the number of colony-forming units (CFU) per gram of stool, expressed in a computer version of scientific notation, where the capital letter “E” indicates the exponent value
- Thus, in the example below, the result for *Bacteroides-Prevotella* group is 4.3E7 (or 4.3 x 10⁷) CFU/g stool, representing 43 million individual organisms per gram of stool



Recall that there is a great deal of variability in terms of specific species present in any one individual gut microbiome, so we should not necessarily expect that a reference population of 100 healthy individuals would have the exactly the same composition of bacteria present. As noted in published research, not all commensal bacteria are present in every individual – or if present, may be present in such small amounts relative to other gut microbial species that the level of sensitivity of the molecular assay cannot detect their presence. For some of the commensal targets in *GI Effects*, a portion of the quintile bar is shaded in grey to represent the proportion of the healthy reference population for whom levels of the specific commensal target were below detection limits. For instance, *Oxalobacter formigenes*, one of the 24 targets on *GI Effects*, has a unique ability to metabolize

oxalates in the gut and therefore colonization with this organism may reduce the risk of calcium oxalate kidney stones with healthy levels associated with a 70% reduced risk of being a recurrent calcium oxalate stone-former. However, we know from the literature that *Oxalobacter* is not normally present in 23% to 54% of healthy adults, which is consistent with our finding that 20% of our healthy cohort had a <detection limit (<DL) result for this bug.



The number of tails (1 tail or 2 tail) for the quintile reporting bar for a specific organism is determined empirically by the underlying distribution of the reference population data (i.e., does the data look like it has 1 tail or 2 tails) and by the clinical associations published in literature (i.e., is there published evidence suggesting that High and/or Low levels of this organism are associated with health or with disease). For example, *Bacteroides vulgatus* could not be detected in 5% of our healthy reference population and therefore the reference range value for this bug is <4.6E9, the 97.5% percentile value. We also know from published literature that high levels of *B. vulgatus* are found to be present in stools of severely autistic children when compared to controls. However, we also know that, compared to healthy controls, low levels of *B. vulgatus* are associated with obesity, T2D, and IBS. Therefore, the quintile reporting bar for *B. vulgatus* is 2-tailed even though we are only able to statistically determine the 97.5% percentile value on the right-hand tail.



Outcomes and General Therapeutic Considerations

Once the profile of the commensal bacteria community for an individual patient is known, interventions can be applied in some cases to modulate it favorably.

Supporting the Patient with Disrupted Commensal Bacteria (Dysbiosis)

Dietary changes:

- Increased intake of fiber and whole, complex carbohydrate and resistant starch⁴¹
 - › Introduce high-fiber foods gradually to avoid exacerbation of GI symptoms⁷⁶

Supplements:

- To promote and sustain beneficial commensals, consider addition of probiotics (live bacteria) and/or prebiotic supplements (products to support indigenous microbiota) to alter commensal bacteria; for example:⁹¹
 - Probiotics
 - › *Lactobacillus*
 - › *Bifidobacteria*
 - Prebiotics
 - › Psyllium
 - › Oat bran
 - › Oligofructose
 - › Xylooligosaccharide
 - › Inulin
 - › Beta-glucan
 - › Arabinogalactan
 - Botanical products may also be used to decrease or modulate gut bacteria

Note that:

- Alterations in the macronutrient composition of the diet, such as in vegan, high-protein, or gluten free diets, have shown changes in the composition of the gut microbiome.^{164,165,166-170}
 - › Short-term dietary changes have been shown to rapidly change the microbiome¹⁷¹
 - Animal-based diets were shown to increase the abundance of bile-tolerant organisms, such as *Bacteroides* (plus other Bacteroidetes genera), with a concomitant relative decrease in saccharolytic bacteria
 - A primarily plant-based diet exhibited higher populations of saccharolytic bacteria, such as *Roseburia* and *F. prausnitzii*, which metabolize dietary plant polysaccharides
- Unbalanced or low commensal bacteria levels have been associated with antibiotic or botanical intake and inappropriate probiotic or prebiotic supplementation



Individual Commensal Bacteria

<i>Bacteroides-Prevotella</i> Group Bacteroidetes phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
	L	L/H	L	L/H		LD
	<ul style="list-style-type: none"> Abundance associated with lower bacterial gene richness in the gut¹⁷² Reduced patterns of <i>Bacteroides</i> reported in IBS and ulcerative colitis;¹⁷³ conversely, other researchers found increased levels in IBD¹⁷⁴⁻¹⁷⁷ When compared with fibromyalgia patients, early RA patients showed less <i>Bacteroides-Prevotella-Porphyrromonas</i>¹⁷⁸ Higher levels associated with excessive weight gain in pregnancy¹⁷⁹ and in obesity^{180,181} Other researchers reported lower <i>Bacteroides</i> (as part of <i>Bacteroides-Prevotella</i> group) in obese subjects compared to lean¹⁸² Ratio of <i>Bacteroides-Prevotella</i> group to other gut bacteria correlated positively and significantly with plasma glucose;¹⁸³ In contrast, some have reported half the <i>Bacteroides</i> abundance in T2DM compared to those with normal glucose tolerance or those with pre-diabetes¹⁸⁴ 					
<i>B. vulgatus</i> Bacteroidetes phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
	L		L/H	L	H	LD
	<ul style="list-style-type: none"> <i>Bacteroides</i> spp. associated with lower bacterial gene richness in the gut¹⁷² Lower levels of <i>B. vulgatus</i> have been seen in IBS patients in comparison to healthy controls¹⁸⁵ Low relative proportions of <i>B. vulgatus</i>, along with high concentrations of <i>Lactobacillus</i> spp. observed in the microbiota of obese children when compared to lean;¹⁸⁶ <i>B. vulgatus</i> also found under-represented in microbiota of type-2 diabetics¹⁸⁷ <i>B. vulgatus</i> found to be present in significantly higher numbers in stools of severely autistic children when compared to controls¹⁸⁸ While increased <i>B. vulgatus</i> prevalence was associated with the genotype of infants at high risk of celiac disease development,¹⁸⁹ another study found that <i>B. vulgatus</i> was more frequently detected in controls than in patients with treated celiac disease ($p < 0.01$)¹⁹⁰ 					
<i>Barnesiella</i> spp. Bacteroidetes phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
			L			
	<ul style="list-style-type: none"> Currently identified as a common but low-abundance genus¹⁹¹ Higher fecal <i>Barnesiella</i> (two logs difference) associated with protective effect against Vancomycin-Resistant Enterococcus in stem cell transplant patients¹⁹² 					
<i>Odoribacter</i> spp. Bacteroidetes phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
		L				
	<ul style="list-style-type: none"> Lower concentrations in humans have been reported in ileal Crohn's and pancolonic ulcerative colitis¹⁹³ 					
<i>Prevotella</i> spp. Bacteroidetes phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
		H	H	H	L	LD
	<ul style="list-style-type: none"> Abundance associated with lower bacterial gene richness in the gut¹⁷² Has been found higher in IBD,^{177,194} in smokers with (and without) IBD (in association with higher <i>Bacteroides</i> and lower <i>F. prausnitzii</i>),¹⁷⁴ and in type-2 diabetes¹⁸⁴ <i>Prevotella</i> found over-represented in new-onset rheumatoid arthritis (RA), when compared to other groups (healthy controls, treated RA patients, and psoriatic arthritis patients, who had higher <i>Bacteroides</i>)¹⁹⁵ Lower levels found in autism¹⁹⁶ and in debilitated aging¹⁹⁷ Family Prevotellaceae exhibited a >6-fold increase in obese subjects when compared to the healthy group, with most of the Prevotellaceae sequences belonging to a single-genus, <i>Prevotella</i>.¹⁹⁸ 					

* LD = Robust levels of this organism associated with Low Diversity of gut bacteria
 HD = Robust levels of this organism associated with High Diversity of gut bacteria

Individual Commensal Bacteria

<i>Anaerotruncus colihominis</i> Firmicutes phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
				L		HD
	<ul style="list-style-type: none"> Abundance associated with higher bacterial gene richness in the gut¹⁷² Inversely associated with BMI and triglycerides^{172,199} 					
<i>Butyrivibrio crossotus</i> Firmicutes phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
				L		HD
	<ul style="list-style-type: none"> Abundance associated with higher bacterial richness in the gut²⁰⁰ Abundance may help protect against weight gain²⁰⁰ 					
<i>Clostridium spp.</i> Firmicutes phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
		L/H			H	LD
	<ul style="list-style-type: none"> Many of its species are associated with lower bacterial gene richness¹⁷² Higher <i>Clostridium</i> counts and increased number of <i>Clostridium</i> species reported in autism,²⁰¹⁻²⁰³ vancomycin (which targets <i>Clostridium</i>) improves symptoms in children with late-onset regressive autism²⁰²⁻²⁰⁴ Both higher¹⁹³ and lower²⁰⁵ abundance of <i>Clostridium</i> has been observed in IBD 					
<i>Coprococcus eutactus</i> Firmicutes phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
	L				L	HD
	<ul style="list-style-type: none"> Abundance associated with greater bacterial gene richness in the gut¹⁷² <i>Coprococcus</i> may be less prevalent in autistic children compared to neurotypical children;¹⁹⁶ may be result of intestinal disaccharidase deficiencies common in autism²⁰⁶ In IBS, reduced abundance reported (in association with elevated <i>Ruminococcus spp.</i>)²⁰⁷ 					
<i>Faecalibacterium prausnitzii</i> Firmicutes phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
		L/H		L		HD
	<ul style="list-style-type: none"> In a healthy gut, represents more than 5% of the total bacterial population²⁰⁸ and is comprised of only one species Abundance associated with higher bacterial gene richness in the gut¹⁷² Controls inflammation through inflammatory-cytokine inhibition,²⁰⁹ lower counts reported in IBD:^{174,193,210} Crohn's disease²¹¹⁻²¹⁵ and ulcerative colitis (UC),²¹⁶ although increases have been noted²¹⁵ Appears to protect against glucose intolerance and type 2 diabetes;¹⁸⁴ possibly due to anti-inflammatory effects¹⁸² and/or positive effects on insulin resistance status^{182,184} 					
<i>Lactobacillus spp.</i> Firmicutes phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
	L/H			H		HD
	<ul style="list-style-type: none"> Abundance associated with higher bacterial gene richness in the gut¹⁷² Studies have reported altered levels in IBS, with some finding higher concentrations²¹⁷⁻²¹⁹ and others finding lower concentrations²²⁰⁻²²² Lower levels reported to correlate with symptom severity in IBS²⁰⁷ Increased levels seen in obese patients compared to lean controls²²³ 					
<i>Pseudoflavonifractor spp.</i> Firmicutes phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
						HD
	<ul style="list-style-type: none"> Abundance associated with higher bacterial gene richness in the gut¹⁷² 					



Individual Commensal Bacteria

Roseburia spp. Firmicutes phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
	L	L	L	L		HD
	<ul style="list-style-type: none"> Abundance associated with higher bacterial gene richness in the gut¹⁹⁹ Less abundant in individuals with IBS, particularly constipation-predominant IBS²²⁴ Counts lower in type 2 diabetics;¹⁸³ trending inversely with plasma glucose¹⁸³ Lower in IBD^{193,216} and early-onset rheumatoid arthritis¹⁷⁸ (as part of decreased <i>E. rectale-C. coccoides</i> group) 					
Ruminococcus spp. Firmicutes phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
	L/H	H			H	LD
	<ul style="list-style-type: none"> Abundance associated with low bacterial gene richness in the gut¹⁷² Human studies have reported that Ruminococcus spp. tend to be more abundant in IBD;²²⁵ active UC,¹⁷⁶ active CD,²²⁶ and ileal CD²²⁷ Levels are variable in IBS, depending on IBS subtype, with some researchers reporting increased concentrations²²⁸ and some finding decreased amounts²²² May be more prevalent in autism²⁰² 					
Veillonella spp. Firmicutes phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
	L/H				L	HD
	<ul style="list-style-type: none"> Abundance associated with higher bacterial gene richness in the gut¹⁷² Imbalances noted in IBS, although findings are mixed: some studies reported higher concentrations in IBS,²¹⁸ in IBS-C,²²¹ IBS-D;²¹⁹ others have reported lower counts²¹⁷ or lower counts weakly correlating with greater symptom severity²⁰⁷ Found less abundant in autistic children compared to neurotypical children¹⁹⁶ 					
Bifidobacterium spp. Actinobacteria phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
	L	L	L	L/H	L	HD
	<ul style="list-style-type: none"> Abundance associated with higher bacterial gene richness in the gut¹⁷² Modulates local and systemic immune responses²²⁹ Abundance lower in IBD^{212,213} Abundance lower in IBS;^{219,230,231} low levels also correlate with symptom severity in IBS²⁰⁷ Lower levels seen in type 2 diabetes,^{187,232,233} pediatric allergy,²³⁴ and autism²³⁵ Increased levels in obese subjects compared to lean/overweight;²³⁶ infants with lower <i>Bifidobacterium</i> may have increased risk for weight gain in childhood²³⁷ Abundance decreases after weight loss²³⁸ and gastric-bypass surgery¹⁸² 					
B. longum Actinobacteria phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
				H		HD
	<ul style="list-style-type: none"> Abundance associated with higher bacterial gene richness in the gut¹⁷² Abundance decreases with weight loss²³⁸ Found increased in obese subjects compared to lean/overweight²³⁶ 					
Collinsella aerofaciens Actinobacteria phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
	L	H		H		
	<ul style="list-style-type: none"> Lower counts reported in IBS;^{222,239} lower levels may correlate with greater severity of IBS symptoms²⁰⁷ Higher concentrations reported in IBD; thought to be result of abnormal host responses to the bacteria¹⁷⁷ <i>Collinsella</i> spp. reported higher in type 2 diabetes;¹⁸⁴ 					

* LD = Robust levels of this organism associated with Low Diversity of gut bacteria
 HD = Robust levels of this organism associated with High Diversity of gut bacteria

Individual Commensal Bacteria

<i>Desulfovibrio piger</i> Proteobacteria phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
	H	H		L	H	
	<ul style="list-style-type: none"> • Reported higher in IBD^{205,240,241} • Sulfate-reducing bacteria higher in constipation-predominant IBS²²⁴ compared with healthy subjects • <i>Desulfovibrio</i> spp. also found higher in autism²⁴² • May be lower in obesity: preschoolers²⁴³ and adults²⁴⁴ 					
<i>Escherichia coli</i> Proteobacteria phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
	H	H		L/H		
	<ul style="list-style-type: none"> • Increased counts reported in inflammatory bowel disease;^{177,205,213,245-247} • Increased levels found in diarrhea-predominant IBS²⁴⁸ • Higher in overweight pregnant women compared to normal weight pregnant women and in women with excessive weight gain during pregnancy²⁴⁹ • Reported to increase with weight loss after gastric bypass, correlating negatively with leptin levels¹⁸² 					
<i>Oxalobacter formigenes</i> Proteobacteria phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
				L		
	<ul style="list-style-type: none"> • Normally present in 46–77% of healthy adults²⁵⁰ • Unique ability to metabolize oxalates in the gut²⁵⁰⁻²⁵² • Dietary oxalate consumption generally increases <i>O. formigenes</i> abundance in controls, but not stone formers²⁵⁰ • Colonization with this bacteria may reduce risk of oxalate stone formation,^{250,252} with healthy levels associated with 70% reduced risk of being recurrent calcium-oxalate stone-former²⁵⁰ 					
<i>Methanobrevibacter smithii</i> Euryarchaeota phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
	L/H	L/H		L/H		HD
	<ul style="list-style-type: none"> • Abundance associated with higher bacterial gene richness in the gut¹⁷² • Lower counts of <i>Methanobrevibacter</i> species reported in human obesity;^{180,253} higher amounts reported in anorexia;²²³ in contrast, one study confirmed a positive association with increased BMI and body fat in methanogen-colonized populations^{254**} • Higher levels linked to IBS-C; reduced levels linked with IBS-D^{255**} • Methanogens found higher in people with colon cancer, colonic polyposis, ulcerative colitis,^{256**} and diverticular disease (sigmoidoscopy enema samples)²⁵⁷ • Some studies have reported lower counts in IBD;^[255,258**] ²⁵⁹ conversely, other have reported increased abundance;^{256**} <p>** Breath-testing study, an accepted indirect measure of gut methanogens; although <i>M. smithii</i> is currently considered the dominant methanogenic archaeon in the gut,^{260,261} other methanogenic bacteria may also be contributors to breath methane.²⁵⁹</p>					
<i>Fusobacterium</i> spp. Fusobacteria phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
		H	H	H		
	<ul style="list-style-type: none"> • Although part of normal human gut flora, species of <i>Fusobacterium</i> strongly associated with numerous diseases, including colorectal cancer (CRC).²⁶²⁻²⁶⁵ appendicitis,²⁶⁶ dental plaque/ periodontal disease,²⁶⁷ hepatic cirrhosis,²⁶⁸ and inflammatory bowel disease^{176,269,270} • <i>Fusobacterium</i> correlates positively with TNF-alpha, suggesting involvement of mucosal inflammation²⁶⁴ • Obese, older subjects with metabolic syndrome demonstrated increased <i>Fusobacterium</i> as compared to younger subjects¹⁹⁹ 					

* LD = Robust levels of this organism associated with Low Diversity of gut bacteria
 HD = Robust levels of this organism associated with High Diversity of gut bacteria



Individual Commensal Bacteria

<i>Akkermansia muciniphila</i> Verrucomicrobia phylum	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity* Association
		L		L	L	HD
	<ul style="list-style-type: none"> • Dominant mucus-layer species;²⁷¹ may represent 3-5% of microbial community in healthy adults²⁷² • Abundance associated with higher bacterial gene richness in the gut¹⁷² • Plays role glucose homeostasis¹⁸⁴ • Abundance inversely correlated with IBD (both Crohn's and UC)^{194,225,273} and appendicitis²⁶⁶ • Abundance inversely correlates with body weight in pregnant women²⁴⁹ and children²⁴³ • Some have reported decreased <i>A. muciniphila</i> in pre-diabetes and decreased Verrocomicrobiae abundance in T2D and pre-diabetes¹⁸⁴ • Lower in autism^{196,335} 					

Firmicutes/Bacteroidetes Ratio

The Biomarker

The two largest phyla making up the gut microbiome in humans are Firmicutes and Bacteroidetes. The relationship of these two large groups, expressed as the Firmicutes/Bacteroidetes ratio, has been associated with a number of pathological conditions.²⁷⁴

- Obesity has been specifically associated with a greater abundance of Firmicutes and/or a drop in Bacteroidetes (i.e., an increase in the ratio); some research, however, has shown no change or even an increase in Bacteroidetes in obesity^{253,274,275}
- Disruptions of metabolic homeostasis²⁷⁶
 - › e.g., type 2 diabetes and non-alcoholic fatty liver disease^{277,278}
- Elevated markers of inflammation such as IL-6²⁷⁹

Biomarker Key Points

The Firmicutes/Bacteroidetes Ratio captures the following members of each class:

Total Firmicutes	Total Bacteroidetes
<i>Anaerotruncus colihominis</i>	<i>Bacteroides-Prevotella</i> Group
<i>Butyrivibrio crossotus</i>	<i>B. vulgatus</i>
<i>Clostridium</i> spp.	<i>Prevotella</i> spp.
<i>Coprococcus eutactus</i>	<i>Barnesiella</i> spp.
<i>Faecalibacterium prausnitzii</i>	<i>Odoribacter</i> spp.
<i>Lactobacillus</i> spp.	
<i>Pseudoflavonifractor</i> spp.	
<i>Roseburia</i> spp.	
<i>Ruminococcus</i> spp.	
<i>Veillonella</i> spp.	

Comparator/Gold Standard Tests

Because of the relative newness of the Firmicutes/Bacteroidetes Ratio as determined by PCR, there are no existing comparator tests.

Interpretation

- The quintile reporting system is used for all commensal bacteria results, including the F/B ratio, based on a clinically characterized healthy reference population.
- The differences observed in the gut microbiota composition of obese and lean subjects in some clinical trials has led to multiple preliminary studies on the use of (beneficial) probiotic bacteria to alter the obese phenotype. These studies have largely shown that probiotic intervention has a beneficial effect, and may lead to novel interventions for overweight or obese human patients.²⁸⁰

Outcomes and General Therapeutic Considerations

Although no specific therapeutic interventions are yet known for a markedly elevated or depressed Firmicutes/Bacteroidetes Ratio, general intervention strategies directed against some of the commensals may impact this ratio and ultimately gut health.

Supporting the Patient with Disrupted F/B Ratio

For imbalances in the F/B Ratio it is best to start by correlating with microbial diversity and abundance, and base treatments on those analytes as noted out-of-range in the four functional pillars.

Bacteriology (Culture)

Traditional bacterial culture complements DNA-based tests such as PCR to provide a more complete survey of a patient's gut microbiota beyond the specific organisms targeted by PCR. Culture methods have established clinical utility and defined parameters that have been long recognized as “gold standard” in traditional clinical diagnostics.

Beneficial Bacteria (Culture)

Cultivable gut bacteria include beneficial bacteria such as *Bifidobacteria*, *Lactobacilli*, and *E. coli*, all of which are thought to exert positive local and systemic effects through their anti-inflammatory and immune-modulating properties. Conversely, impaired commensal bacterial populations have been associated with intestinal and chronic diseases. *Bifidobacteria*, *Lactobacilli*, and *E. coli* are also targeted by PCR to provide more specific quantitation of their levels.

Biomarker Key Points

- *Lactobacillus* species and *E. coli* are facultative anaerobes, capable of surviving in an environment with limited oxygen.
- *Bifidobacteria* are obligate anaerobes, and must be grown in anaerobic chambers.

Culture is required for determining therapeutic interventions such as sensitivities to pharmaceutical or botanical antibiotics.

Patient Populations of Interest

Bacterial culture remains of great utility in evaluating patients with symptoms of possible gastrointestinal infection or dysbiosis.





Comparator/Gold Standard Tests

Bacterial culture is the gold standard for identification of populations of cultivable organisms.


Interpretation

Results of bacteriology culture are reported as “No Growth,” or growth in one of three categories of bacteria, using a color-coding system:

- No growth (white)
- Non-pathogen (green)
- Potential pathogens (yellow)
- Known pathogens (red)

Microbiology Legend			
NG	NP	PP	P
			
No Growth	Non-Pathogen	Potential Pathogen	Pathogen

The quintile reporting system is not used for this test; instead, growth is reported as present at one of four semi-quantitative levels of abundance.



Outcomes and General Therapeutic Considerations

Patients with low growth of beneficial bacteria may benefit from dietary manipulations and supplements similar to those discussed in the section on commensal bacteria, including pre- and probiotics and elimination of antibiotics when feasible.

Additional Bacteria (Culture)

Bacterial culture of fecal material may also yield additional organisms of interest in the assessment of the gut microbiome. Many such organisms constitute part of the normal aerobic or commensal flora, and are not recognized as major pathogens. They are typically readily cultured for identification when they occur at clinically significant levels.

Additional bacteria also include organisms that may be potential pathogens (PP on the report form), or “opportunistic” pathogens. Not usually pathogenic, these organisms may overgrow during periods of perturbation in the gut environment and their presence on culture may be an indicator of imbalance or dysbiosis. The significance of these organisms must be interpreted in the context of a specific patient's clinical presentation (e.g. symptoms, immunosuppression, etc.).

Biomarker Key Points

Additional or opportunistic pathogens are typically restrained and controlled by balanced levels of commensal organisms, but their overgrowth may occur when commensal bacterial populations are impaired by:

- Infection with overt pathogens or parasites
- Poor diet
- Antibiotic use
- Lowered gut immunity

Patient Populations of Interest

Patients whose management might benefit from testing for Additional Bacteria include those with symptoms of unexplained or persistent diarrhea, especially those with known or suspected imbalance of the normal gut flora (dysbiosis), and those who are immunosuppressed.

Comparator/Gold Standard Tests

This test uses bacterial culture. Identification of cultured bacteria is via Vitek-MS using Matrix Assisted Laser Desorption Ionization Time-Of-Flight (MALDI-TOF). The MALDI-TOF mass spectrometry platform utilized for the rapid identification of bacteria and yeast from pure cultures on the *GI Effects Comprehensive Profile* report relies on the most extensive FDA-cleared library of microbial targets available on the market, which can accurately identify approximately 200 different additional bacterial species.

Interpretation

Additional bacteria are reported using the same conventions as for Bacteriology Culture, above.

- No growth (white)
- Non-pathogen (green)
- Potential pathogens (yellow)
- Known pathogens (red)

Microbiology Legend			
NG	NP	PP	P
No Growth	Non-Pathogen	Potential Pathogen	Pathogen

The quintile reporting system is not used for this test; instead, growth of any additional bacteria is reported as present at one of four semi-quantitative levels of abundance, as shown in the following example:

Additional Bacteria	+3	NP	
<i>Alpha</i> haemolytic streptococcus	+3	NP	
<i>Gamma</i> haemolytic streptococcus	+4	PP	
<i>Citrobacter freundii</i>	+2	NP	

The presence of Potentially Pathogenic organisms (PP) at higher levels, or of Pathogenic bacteria (P) at any level, should trigger increased concern. Such levels may be of clinical relevance in patients with bacterial gastroenteritis including:

- Travelers' diarrhea²⁸¹⁻²⁸⁵
- Food poisoning²⁸⁶
- IBD²⁸⁷
- IBS²⁸⁸

Presence of opportunistic or potentially pathogenic organisms may also indicate intestinal microflora imbalance, or dysbiosis, poor diet, antibiotic use, or lowered gut immunity.

Outcomes and General Therapeutic Considerations

Once identified, pathogenic and potentially pathogenic/opportunistic bacteria may be decreased through pharmacological and botanical treatments, and by restoring a healthy intestinal microflora balance.

Pharmacological treatments are best utilized subsequent to standard methods of diagnosis, including sensitivity profiles. When opportunistic bacteria are found, they are automatically cultured and the sensitivity of pharmaceuticals and active botanical ingredients are assessed. Agents marked Sensitive or "S" have been shown to be effective treatments in culture/sensitivity studies and may be used as indicated.



Mycology (Culture) The Biomarker

The fungal kingdom includes yeasts and molds. Yeasts are single-celled organisms in the fungal kingdom; 1500 species have been described to date. All humans house fungal colonies in their colons, typically without detriment. Yeasts are likely to be normal members of the human microflora and certain strains (e.g., *Saccharomyces boulardii*) have been shown to:

- Reduce symptoms of diarrhea in children
- Prevent reinfection with *C. difficile*
- Reduce bowel movements in patients with diarrhea-predominant IBS
- Reduce the incidence of antibiotic-, travelers'- or HIV-associated diarrhea

Biomarker Key Points

- Fungi are detected in fecal samples by standard culture techniques
- Pathogenic and potentially pathogenic (opportunistic) fungi are associated with gastrointestinal symptoms, especially in immune-compromised people
- Fungal infections can produce imbalances of GI microorganisms (dysbiosis)

Patient Populations of Interest

Patients at increased risk of significant fungal infections (mycosis) include people with:

- Compromised immunity
 - › Those on corticosteroids
 - › Those with diabetes
 - › The very young or very old
 - › Yeast overgrowth syndrome (clinical presentation does not have the severity described in conventionally recognized fungal infections)

Comparator/Gold Standard Tests


This test uses fungal culture, which is the gold standard. Identification of cultured yeast/fungi is via Vitek-MS (MALDI-TOF). The MALDI-TOF mass spectrometry platform utilized for the rapid identification of bacteria and yeast from pure cultures on the *GI Effects Comprehensive Profile* report relies on the most extensive FDA-cleared library of microbial targets available on the market.

Interpretation

Fungal culture is reported as Mycology (Culture) on the test report, using the same conventions as for Bacteriology Culture, above.

- No growth (white)
- Non-pathogen (green)
- Potential pathogens (yellow)
- Known pathogens (red)

The quintile reporting system is not used for this test; instead, growth is reported as present at one of four semi-quantitative levels of abundance. The example shows growth of a potential pathogen (PP) in the second quadrant of growth.



When fungi are found, they are automatically cultured and the sensitivities of pharmaceuticals and active botanical ingredients are assessed. Agents marked Sensitive or “S” have been shown to be effective treatments in culture/sensitivity studies and may be used as indicated.

Outcomes and General Therapeutic Considerations

Once identified, treatment for pathological yeast species and yeast overgrowth should be directed at eliminating the offending species and restoring microbial balance. Treatments may be pharmaceutical or botanical. Pharmaceutical treatments should be used following standard methods of diagnosis.

Supporting the Patient with Abnormal Fungal Culture

In supporting patients with yeast problems, consider the degree of infection, the patient’s overall immune status, diet (carbohydrate/sugar intake), and offer support for a potentially impaired immune system.

Parasitology The Biomarker

Though prevalence data for intestinal parasitic infection in the US is limited, one survey found one-third of 5792 fecal specimens to be positive for parasites, with positivity peaking seasonally between July and October.²⁸⁹ According to the American Association for Clinical Chemistry (AACC), the most common parasites in the United States include: *Cryptosporidium*, *Entamoeba histolytica*, and *Giardia*.

A single fecal specimen tested for parasites by O&P may detect approximately 90% of GI parasite infections.⁸⁹

The parasitology test uses two complementary methodologies:

1. Microscopic examination of fecal specimens for ova and parasites (O&P), the gold standard of diagnosis for many parasites
2. Enzyme immunoassay (EIA) for the identification of *Cryptosporidium*, *Entamoeba histolytica*, and *Giardia lamblia*

Biomarker Key Points

- EIA is a biochemistry-based test that detects immunogenic macromolecules such as toxins or organism-specific antigens
- EIA is widely recognized for its diagnostic utility for detection of pathogenic antigens

Patient Populations of Interest

Screening for parasites is appropriate in patients presenting with:

- Excessive or persistent diarrhea
- Stools containing blood or mucous
- Severe abdominal pain (chronic or subacute)
- Nausea and vomiting

Because of the diversity of presentations in patients with GI parasites, it may be useful to test for parasites in patients with such symptoms who:

- Have unexplained persistent headache and fatigue
- Have been exposed to a parasitic outbreak at daycare or school
- Have traveled outside of the U.S.
- Have consumed untreated water

Comparator/Gold Standard Tests

Microscopic examination of fecal samples for ova and parasites is the gold standard test for such examinations, and continues to have the highest proven diagnostic and clinical utility for parasite detection.

Although the examination of at least three samples on at least 3 separate days remains the recommendation of the Centers for Disease Control, some literature suggests that approximately 90% of enteric parasite infections may be detected in a single stool sample collected for O&P examination, with small increases in sensitivity and negative predictive values for additional samples.⁸⁹ For patients where parasitic infection is suspected or needs to be excluded, at least three samples on three separate days should be submitted for evaluation.

Interpretation

The normal result for O&P on microscopic exam is “negative” or “no organisms seen.” Parasitology EIA tests are reported as “Negative” (in range) or “Positive” (out of range) on the report, with each EIA-evaluated organism listed separately (*Cryptosporidium*, *E. histolytica*, and *G. lamblia*).

Outcomes and General Therapeutic Considerations

Because parasitic infections typically require pharmaceutical treatment with anti-parasite medications, an accurate and timely diagnosis is essential. In a patient with a high degree of clinical suspicion for parasites, a positive test warrants appropriate interventions.

The Table below lists common fecal parasitic pathogens along with symptoms and therapeutic considerations.



Common Parasitic Protozoans

Parasitic Protozoans	Symptoms	Therapeutic Considerations
<p><i>Blastocystis hominis</i> Transmission: Unknown; consider contaminated food or water, or exposure to infected animals suspected.</p>	<p>Watery or loose stools, diarrhea, abdominal pain, anal itching, weight loss, constipation, fatigue, and excess flatulence have been reported in persons with <i>Blastocystis</i> infection. Many people are asymptomatic.</p>	<p>The clinical significance of <i>Blastocystis</i> spp. is controversial, although there is increasing evidence that it may be a pathogen in some individuals with symptoms meeting criteria for Irritable Bowel Syndrome.⁹⁰ For additional information, see www.dpd.cdc.gov/dpdx/HTML/Blastocystis.htm</p>
<p><i>Cryptosporidium</i> spp. Transmission: Fecal contamination of food or water, including swimming pools and municipal water supplies.</p>	<p>Patients will be symptomatic or present with diarrhea varying from mild to severe, abdominal cramping, weight loss, anorexia, nausea, vomiting, flatulence, malaise, and mild fever.</p>	<p>Most people who have healthy immune systems will recover without treatment. Diarrhea can be managed by drinking plenty of fluids to prevent dehydration. Immunosuppression increases infection severity. See www.dpd.cdc.gov/dpdx/HTML/Cryptosporidiosis.htm</p>
<p><i>Dientamoeba fragilis</i> Transmission: Unknown; often associated with pinworm infection or fecal contamination.</p>	<p>Patients will be asymptomatic or present with diarrhea, nausea, and vomiting; abdominal tenderness is possible.</p>	<p>There is no consensus as to best clinical practice; goal is eradication of parasite. See www.dpd.cdc.gov/dpdx/HTML/Dientamoeba.htm</p>
<p><i>Entamoeba coli</i>, <i>Entamoeba histolytica</i> and <i>Entamoeba dispar</i> Transmission: Contaminated food or water, pets, sexual contact.</p> <p>More common in people who live in tropical areas with poor sanitary conditions. It can be a pathogenic amoeba.</p>	<p>Several protozoan species in the genus <i>Entamoeba</i> colonize humans, but not all of them are associated with disease. <i>Entamoeba histolytica</i> is well recognized as a pathogenic amoeba, associated with intestinal and extra-intestinal infections. Only about 10% to 20% of people who are infected with <i>E. histolytica</i> become sick. A severe form of <i>E. histolytica</i> is associated with stomach pain, bloody stools, and fever (may resemble ulcerative colitis). <i>E. dispar</i> is non-pathogenic.</p>	<p>Only one antibiotic is used in non-symptomatic <i>E. histolytica</i> infection; two antibiotics if patients are symptomatic. See www.dpd.cdc.gov/dpdx/HTML/IntestinalAmebae.htm and www.dpd.cdc.gov/dpdx/HTML/Amebiasis.htm</p>
<p><i>Giardia lamblia</i> Transmission: Contaminated water, food, or fecal-oral transmission. <i>G. lamblia</i> is the leading cause of intestinal parasitic infection in the US. See www.dpd.cdc.gov/dpdx/HTML/Giardiasis.htm</p>	<p>Patients can be asymptomatic. If symptomatic, will present with acute to chronic diarrhea with bloating, intestinal malabsorption, and steatorrhea. Giardiasis has been associated with agammaglobulinemia, chronic pancreatitis, achlorhydria, and cystic fibrosis.</p>	<p>Several prescription drugs are available to treat giardiasis. See emedicine.medscape.com/article/176400-clinical</p>



Bacteria Sensitivity

When bacterial culture yields pathogenic or potentially pathogenic organisms, the Bacteria Sensitivity section reports the results of in vitro testing for susceptibility.

- **Prescriptive Agents** lists detected organisms along with their relative sensitivity to prescription antibiotics or antimicrobials. Conventional definitions of “sensitive,” “intermediate,” and “resistant” organisms are used.
- **Natural Agents** lists detected organisms along with their relative degree of growth inhibition by herbal and other natural substances. “High” inhibition indicates a greater ability of the substance to limit microbial growth, while “Low” inhibition suggests less ability to limit growth.

Mycology Sensitivity

When fungal culture of stool yields fungal organisms, the Mycology Sensitivity section reports the results of in vitro testing for susceptibility. Prescriptive and Natural Agents are listed and categorized similarly to those reported for Bacteria Sensitivity.

- **Prescriptive Agents** lists detected organisms along with their relative sensitivity to prescription antibiotics or antimicrobials. Conventional definitions of “sensitive,” “intermediate,” and “resistant” organisms are used.
- **Natural Agents** lists detected organisms along with their relative degree of growth inhibition by herbal and other natural substances. “High” inhibition indicates a greater ability of the substance to limit microbial growth, while “Low” inhibition suggests less ability to limit growth.

Additional Tests

Several additional tests have long been used in the analysis of stool as a diagnostic analyte. These include the color, the consistency, and the presence or absence of occult blood.

Color

Stool color is primarily associated with diet and medication use, though it may also be an indicator of various GI health conditions.

Consistency

Stool consistency may vary from significantly hard to watery. The technical ability to recover diagnostic biomarkers from stool may be influenced by extremes of consistency. Consistency is self-reported by the patient.

Occult Blood

The term “occult” in this context simply means blood not evident to the naked eye, that is, blood present in microscopic quantities only. Normally, stools should be entirely free of blood.

The Hemosure diagnostic kit uses fecal immunochemical testing (FIT); it has higher specificity than the common guaiac test because of its use of mono- and polyclonal antibodies specific to human hemoglobin.

FIT-based diagnostics have been recommended by the American College of Gastroenterology as the preferred test for colorectal cancer screening/detection. Further investigation is warranted in the presence of a positive test.

There are no drug or dietary restrictions prior to collecting the sample in relationship to FIT; vitamins or foods do not affect the FIT.

Pathogenic Bacteria (Add-on) The Biomarker

Pathogenic bacteria are those organisms known to cause distinct human disease processes, and in that way they differ from normal or opportunistic organisms in the human GI microbiome.

Pathogenic bacteria are detected on this test by enzyme immunoassay (EIA), which is the standard approach to diagnosis of pathogens.

Biomarker Key Points

Organisms for which EIA testing is done are:

- *Helicobacter pylori* (*H. pylori*)
- *Campylobacter* species
- Shiga Toxin-producing *E. coli*
- *Clostridium difficile* (*C. difficile*)

Patient Populations of Interest

Pathogenic Bacteria EIA testing will achieve its optimum utility when used in the context of an appropriate differential diagnosis that considers patient symptoms and produces a high index of suspicion for at least one clinically known syndrome or symptom complex. Testing of non-symptomatic patients is not recommended.

For each organism, symptomatology is as follows:

- *H. pylori*
 - › Although prevalence of the organism is 35 to 40% among US adults, most patients remain asymptomatic²⁹⁰ (“colonized” and not “infected”) and therefore do not require testing to document the organism’s presence
 - › Indications for testing include upper GI symptoms or pathology such as:
 - Gastritis
 - Duodenal and peptic ulcer disease
 - Gastric lymphoma
 - Gastric cancer (in the patient or a relative)
- *Campylobacter* species
 - › Asymptomatic carriage is common; testing of asymptomatic patients may be indicated during outbreaks

- › Indications for testing include:
 - Diarrhea (often bloody)
 - Abdominal cramps
 - Fever
- › Clinically apparent illness typically lasts about one week, but may be serious in immune-compromised hosts
- › An earlier culture-independent detection method using a DNA probe was less specific, detecting any *Campylobacter* species, not necessarily pathogenic ones
- Shiga toxin-producing *E. coli*
 - › Asymptomatic carriage can occur
 - › When present, symptoms are variable and may include:
 - Abdominal cramping
 - Water or bloody diarrhea
 - Vomiting
 - › EIA confirms the presence of the pathogenic shiga toxin
 - › Instead of routine screening, recommendations for testing now concentrate on clinically relevant populations; specifically, patients with significant diarrhea
- *C. difficile*
 - › Asymptomatic colonization can occur; true infection is defined as presence of the organism and/or its toxin in the context of a symptomatic patient
 - › When present, symptoms include:
 - Cramping
 - Lower abdominal pain/tenderness
 - Fever
 - Watery diarrhea
 - Nausea
 - Loss of appetite
- *C. difficile* is a well-known cause of antibiotic-associated diarrhea

Comparator/Gold Standard Tests

There is substantial peer reviewed literature indicating the diagnostic utility of the detection of a pathogenic antigen by EIA.

- EIA provides actionable clinical information

Interpretation

The positive finding for any pathogen on the report is considered significant.

Outcomes and General Therapeutic Considerations

Therapeutic intervention is warranted in any patient in whom the practitioner has a high clinical index of suspicion and in whom a diagnosis of any of these four pathogens is made.

The table below summarizes the organisms, their typical symptoms, and appropriate therapeutic considerations.



Pathogenic Bacteria by EIA

Pathogen	Symptoms	Therapeutic Considerations
<p><i>Campylobacter</i> spp. Transmission: Consumption of contaminated food (particularly poultry), water, or contact with infected animals (particularly cats and puppies).</p>	<p>Symptoms include: diarrhea (often bloody), abdominal cramps, and fever. Illness typically lasts one week and may be serious in immunocompromised patients. <i>Campylobacter jejuni</i> infection has been associated with the onset of Guillain-Barré syndrome (GBS).</p>	<p>Patients generally recover without any specific treatment. They should drink extra fluids throughout the duration of diarrheal episodes. CDC review is at www.cdc.gov/nczved/divisions/dfbmd/diseases/campylobacter/</p>
<p><i>Clostridium difficile</i> Transmission: Shed in feces.</p>	<p>Symptoms include: cramping, lower abdominal pain/tenderness, fever, watery diarrhea, loss of appetite, and nausea. It is a common cause of antibiotic-associated diarrhea (AAD).</p>	<p>In about 20% of patients, <i>C. diff</i> infection will resolve within 2 to 3 days of discontinuing the antibiotic to which the patient was previously exposed. CDC review is at www.cdc.gov/HAI/organisms/cdiff/Cdiff_faqs_HCP.html#a9</p>
<p>Shiga toxin <i>E. coli</i> Transmission: Consumption of contaminated food, unpasteurized (raw) milk, water that has not been disinfected, contact with cattle or the feces of infected people. STEC stands for shiga toxin-producing <i>Escherchia coli</i> (STEC).</p>	<p>The symptoms vary but may include abdominal cramping, watery or bloody diarrhea, and vomiting.</p>	<p>The infection can be self-limiting. Rehydrate, and consider pre-and probiotics to support infection resolution. Treatment is based on the site and severity of infection, and STEC status. CDC review at www.cdc.gov/ecoli/index.html. Also see www.cdc.gov/ecoli/clinicians.html</p>
<p><i>Helicobacter pylori</i> Transmission: Is incompletely characterized. Person-to-person transmission is most commonly implicated with fecal/oral, oral/oral, or gastric/oral pathways.</p>	<p>Symptoms include: acute gastritis with abdominal pain, nausea, and vomiting. Non-ulcer dyspepsia is common. Development of severe <i>H. pylori</i> disease is partially determined by the virulence of the infecting strain.</p>	<p>Conventional recommendation is polypharmacy: antibiotics and proton-pump inhibitors (PPIs). See www.acg.gi.org. Botanical anti-<i>H. pylori</i> formulas may prove helpful. See www.cdc.gov/ulcer/keytocure.htm.</p>

References

1. Flint HJ, Scott KP, Louis P, Duncan SH. The role of the gut microbiota in nutrition and health. *Nat Rev Gastroenterol Hepatol*. Oct 2012;9(10):577-589.
2. Nicholson JK, Holmes E, Wilson ID. Gut microorganisms, mammalian metabolism and personalized health care. *Nature reviews. Microbiology*. May 2005;3(5):431-438.
3. Plotnikoff GA, Riley D. The human microbiome. *Global advances in health and medicine : improving healthcare outcomes worldwide*. May 2014;3(3):4-5.
4. Structure, function and diversity of the healthy human microbiome. *Nature*. Jun 14 2012;486(7402):207-214.
5. Goepf J FE, McBride T, Landis D. Frequency of Abnormal Fecal Biomarkers in Irritable Bowel Syndrome. *Global Advances in Health and Medicine*. 2014;3(3):9-15.
6. Grigsby B, Rodriguez-Rilo H, Khan K. Antioxidants and chronic pancreatitis: theory of oxidative stress and trials of antioxidant therapy. *Dig Dis Sci*. Apr 2012;57(4):835-841.
7. Farrell JJ. Overview and diagnosis of malabsorption syndrome. *Seminars in gastrointestinal disease*. Oct 2002;13(4):182-190.
8. Loser C, Mollgaard A, Folsch UR. Faecal elastase 1: a novel, highly sensitive, and specific tubeless pancreatic function test. *Gut*. 10/1996 1996;39(4):580-586.
9. Cohen JR, Schall JI, Ittenbach RF, Zemel BS, Stallings VA. Fecal elastase: pancreatic status verification and influence on nutritional status in children with cystic fibrosis. *Journal of pediatric gastroenterology and nutrition*. Apr 2005;40(4):438-444.
10. Apte MV, Pirola RC, Wilson JS. Mechanisms of alcoholic pancreatitis. *J Gastroenterol Hepatol*. Dec 2010;25(12):1816-1826.
11. Burkovskaia VA, Beloborodova EI, Glinskaia ON, et al. [Exocrine function of the pancreas in patients with ulcerative colitis]. *Klinicheskaia meditsina*. 2010;88(1):46-49.
12. Evans KE, Leeds JS, Morley S, Sanders DS. Pancreatic insufficiency in adult celiac disease: do patients require long-term enzyme supplementation? *Digestive diseases and sciences*. Oct 2010;55(10):2999-3004.
13. Hardt PD, Bretz L, Krauss A, et al. Pathological pancreatic exocrine function and duct morphology in patients with cholelithiasis. *Digestive diseases and sciences*. Mar 2001;46(3):536-539.
14. Herzig KH, Purhonen AK, Rasanen KM, et al. Fecal pancreatic elastase-1 levels in older individuals without known gastrointestinal diseases or diabetes mellitus. *BMC geriatrics*. 2011;11:4.
15. Chowdhury RS, Forsmark CE. Review article: Pancreatic function testing. *Aliment.Pharmacol. Ther*. 3/15/2003 2003;17(6):733-750.
16. Ayling RM. New faecal tests in gastroenterology. *Annals of clinical biochemistry*. Jan 2012;49(Pt 1):44-54.
17. Muñoz JED. Diagnosis of chronic pancreatitis: Functional testing. *Best Practice & Research Clinical Gastroenterology*. Vol 242010:233-241.
18. Wali PD, Loveridge-Lenza B, He Z, Horvath K. Comparison of fecal elastase-1 and pancreatic function testing in children. *Journal of pediatric gastroenterology and nutrition*. Feb 2012;54(2):277-280.
19. Walkowiak J, Madry E, Lisowska A, et al. Adaptive changes of pancreatic protease secretion to a short-term vegan diet: influence of reduced intake and modification of protein. *The British journal of nutrition*. Jan 2012;107(2):272-276.
20. Walkowiak J, Wadolowska L, Szafarska-Poplawska A, Lisowska A, Bugajewska A, Przyslawski J. The elimination of meat from the diet selectively decreases pancreatic elastase secretion. *The British journal of nutrition*. Jul 2007;98(1):154-158.
21. Weitgasser R, Abrahamian H, Clodi M, Fortunat W, Hammer H. [Position paper: exocrine pancreatic insufficiency and diabetes mellitus.]. *Wien Klin Wochenschr*. Dec 19 2012.
22. Dominguez-Munoz JE. Pancreatic enzyme replacement therapy for pancreatic exocrine insufficiency: when is it indicated, what is the goal and how to do it? *Adv Med Sci*. 2011;56(1):1-5.
23. Drossman DA, Dumitrascu DL. Rome III: New standard for functional gastrointestinal disorders. *J Gastrointestin Liver Dis*. Sep 2006;15(3):237-241.
24. Hammer HF. Pancreatic exocrine insufficiency: diagnostic evaluation and replacement therapy with pancreatic enzymes. *Digestive diseases (Basel, Switzerland)*. 2010;28(2):339-343.
25. Dominguez-Munoz JE. Pancreatic enzyme therapy for pancreatic exocrine insufficiency. *Gastroenterol Hepatol (N Y)*. Jun 2011;7(6):401-403.
26. Fieker A, Philpott J, Armand M. Enzyme replacement therapy for pancreatic insufficiency: present and future. *Clinical and experimental gastroenterology*. 2011;4:55-73.
27. Evans KE, Leeds JS, Morley S, Sanders DS. Pancreatic insufficiency in adult celiac disease: do patients require long-term enzyme supplementation? *Dig.Dis.Sci*. 10/2010 2010;55(10):2999-3004.
28. Nyangale EP, Mottram DS, Gibson GR. Gut microbial activity, implications for health and disease: the potential role of metabolite analysis. *Journal of proteome research*. Dec 7 2012;11(12):5573-5585.
29. Windey K, De Preter V, Verbeke K. Relevance of protein fermentation to gut health. *Molecular nutrition & food research*. Jan 2012;56(1):184-196.
30. Macfarlane GT, Macfarlane S. Bacteria, colonic fermentation, and gastrointestinal health. *J AOAC Int*. Jan-Feb 2012;95(1):50-60.
31. Hughes R, Magee EA, Bingham S. Protein degradation in the large intestine: relevance to colorectal cancer. *Current issues in intestinal microbiology*. Sep 2000;1(2):51-58.
32. Frenhani PB, Burini RC. [Mechanism of action and control in the digestion of proteins and peptides in humans]. *Arq Gastroenterol*. Jul-Sep 1999;36(3):139-147.
33. Champagne ET. Low gastric hydrochloric acid secretion and mineral bioavailability. *Advances in experimental medicine and biology*. 1989;249:173-184.
34. Pali-Schöll I, Herzog R, Wallmann J, et al. Antacids and dietary supplements with an influence on the gastric pH increase the risk for food sensitization. *Clin Exp Allergy*. Jul 2010;40(7):1091-1098.
35. Untersmayr E, Jensen-Jarolim E. The role of protein digestibility and antacids on food allergy outcomes. *The Journal of allergy and clinical immunology*. Jun 2008;121(6):1301-1308; quiz 1309-1310.
36. Pali-Schöll I, Jensen-Jarolim E. . Anti-acid medication as a risk factor for food allergy. *Allergy*. 2011;66(4):469-477.
37. Nakamura T, Tabeke K, Terada A, et al. Short-chain carboxylic acid in the feces in patients with pancreatic insufficiency. *Acta Gastroenterol Belg*. Sep-Dec 1993;56(5-6):326-331.
38. Zarling EJ, Ruchim MA. Protein origin of the volatile fatty acids isobutyrate and isovalerate in human stool. *J Lab Clin Med*. May 1987;109(5):566-570.
39. Russell WR, Gratz SW, Duncan SH, et al. High-protein, reduced-carbohydrate weight-loss diets promote metabolite profiles likely to be detrimental to colonic health. *The American journal of clinical nutrition*. May 2011;93(5):1062-1072.
40. Rasmussen HS, Holtug K, Mortensen PB. Degradation of amino acids to short-chain fatty acids in humans. An in vitro study. *Scandinavian journal of gastroenterology*. Mar 1988;23(2):178-182.



41. Ruiz A. Overview of Malabsorption. *Merck Manual Online* 2014; http://www.merckmanuals.com/professional/gastrointestinal_disorders/malabsorption_syndromes/overview_of_malabsorption.html?qt=fecal%20nitrogen&alt=sh. Accessed June 17, 2014.
42. Swanson KS, Grieshop CM, Flickinger EA, et al. Fructooligosaccharides and *Lactobacillus acidophilus* modify bowel function and protein catabolites excreted by healthy humans. *The Journal of nutrition*. Oct 2002;132(10):3042-3050.
43. Lee MJ, Crook T, Noel C, Levinson UM. Detergent extraction and enzymatic analysis for fecal long-chain fatty acids, triglycerides, and cholesterol. *Clin Chem*. Dec 1994;40(12):2230-2234.
44. De Smet E, Mensink RP, Plat J. Effects of plant sterols and stanols on intestinal cholesterol metabolism: suggested mechanisms from past to present. *Molecular nutrition & food research*. Jul 2012;56(7):1058-1072.
45. Ahnen DJ, Guerciolini R, Hauptman J, Blotner S, Woods CJ, Wargovich MJ. Effect of orlistat on fecal fat, fecal biliary acids, and colonic cell proliferation in obese subjects. *Clinical gastroenterology and hepatology: the official clinical practice journal of the American Gastroenterological Association*. Nov 2007;5(11):1291-1299.
46. Miettinen TA, Vuoristo M, Nissinen M, Jarvinen HJ, Gylling H. Serum, biliary, and fecal cholesterol and plant sterols in colectomized patients before and during consumption of stanol ester margarine. *Am J Clin Nutr*. May 2000;71(5):1095-1102.
47. Zhi J, Melia AT, Guerciolini R, et al. Retrospective population-based analysis of the dose-response (fecal fat excretion) relationship of orlistat in normal and obese volunteers. *Clinical pharmacology and therapeutics*. Jul 1994;56(1):82-85.
48. Krums LM, Parfenov AI, Ekisenina NI. [Disorders of lipid metabolism in patients with chronic diseases of the small intestine]. *Klinicheskaia meditsina*. Nov 1990;68(11):54-57.
49. Cerda JJ, Artnak EJ. Nutritional aspects of malabsorption syndromes. *Comprehensive therapy*. Nov 1983;9(11):35-46.
50. Tormo R. Increased fecal fat and protein intolerance. *Journal of pediatric gastroenterology and nutrition*. Apr 2014;58(4):e42-43.
51. Brody T. Nutritional Biochemistry. 2nd. ed. *San Diego: Academic Press*. 1999:1006.
52. Belli DC, Levy E, Darling P, et al. Taurine improves the absorption of a fat meal in patients with cystic fibrosis. *Pediatrics*. Oct 1987;80(4):517-523.
53. Carrasco S, Codoceo R, Prieto G, Lama R, Polanco I. Effect of taurine supplements on growth, fat absorption and bile acid on cystic fibrosis. *Acta Universitatis Carolinae. Medica*. 1990;36(1-4):152-156.
54. Gaskin KJ. Nutritional care in children with cystic fibrosis: are our patients becoming better? *European journal of clinical nutrition*. May 2013;67(5):558-564.
55. Bunn SK, Bisset WM, Main MJ, Golden BE. Fecal calprotectin as a measure of disease activity in childhood inflammatory bowel disease. *J Pediatr Gastroenterol Nutr*. 2001 2001;32(2):171-177.
56. Roseth AG, Aadland E, Jahnsen J, Raknerud N. Assessment of disease activity in ulcerative colitis by faecal calprotectin, a novel granulocyte marker protein. *Digestion*. 1997 1997;58(2):176-180.
57. United States Food and Drug Administration. PhiCal™ Fecal Calprotectin Immunoassay - K050007. *Medical Devices Approvals and Clearances* 2013; <http://www.fda.gov/MedicalDevices/ProductsandMedicalProcedures/DeviceApprovalsandClearances/Recently-ApprovedDevices/ucm078192.htm>. Accessed June 17, 2014.
58. Tibble JA, Bjarnason I. Non-invasive investigation of inflammatory bowel disease. *World J Gastroenterol*. Aug 2001;7(4):460-465.
59. Tibble JA, Sigthorsson G, Foster R, Forgacs I, Bjarnason I. Use of surrogate markers of inflammation and Rome criteria to distinguish organic from nonorganic intestinal disease. *Gastroenterology*. Aug 2002;123(2):450-460.
60. Guidelines—Rome III Diagnostic Criteria for Functional Gastrointestinal Disorders. *J Gastrointestin Liver Dis*. Sep 2006;15(3):307-312.
61. Maiden L, Takeuchi K, Baur R, et al. Selective white cell apheresis reduces relapse rates in patients with IBD at significant risk of clinical relapse. *Inflammatory bowel diseases*. Oct 2008;14(10):1413-1418.
62. Fagerhol MK. Calprotectin, a faecal marker of organic gastrointestinal abnormality. *Lancet*. Nov 25 2000;356(9244):1783-1784.
63. Tibble JA, Teahon K, Thjodleifsson B, et al. A simple method for assessing intestinal inflammation in Crohn's disease. *Gut*. 2000 2000;47(4):506-513.
64. Tibble JA, Sigthorsson G, Bridger S, Fagerhol MK, Bjarnason I. Surrogate markers of intestinal inflammation are predictive of relapse in patients with inflammatory bowel disease. *Gastroenterology*. Jul 2000;119(1):15-22.
65. Bjarnason I, Sherwood R. Fecal calprotectin: a significant step in the noninvasive assessment of intestinal inflammation. *Journal of pediatric gastroenterology and nutrition*. 2001 2001;33(1):11-13.
66. Rothenberg ME, Mishra A, Brandt EB, Hogan SP. Gastrointestinal eosinophils. *Immunol Rev*. 2001 2001;179:139-155.
67. Carlson M, Raab Y, Peterson C, Hallgren R, Venge P. Increased intraluminal release of eosinophil granule proteins EPO, ECP, EPX, and cytokines in ulcerative colitis and proctitis in segmental perfusion. *Am J Gastroenterol*. 7/1999 1999;94(7):1876-1883.
68. Peterson CG, Eklund E, Taha Y, Raab Y, Carlson M. A new method for the quantification of neutrophil and eosinophil cationic proteins in feces: establishment of normal levels and clinical application in patients with inflammatory bowel disease. *Am J Gastroenterol*. 2002 2002;97(7):1755-1762.
69. Majamaa H, Laine S, Miettinen A. Eosinophil protein X and eosinophil cationic protein as indicators of intestinal inflammation in infants with atopic eczema and food allergy. *Clin Exp Allergy*. 11/1999 1999;29(11):1502-1506.
70. Levy AM, Gleich GJ, Sandborn WJ, Tremaine WJ, Steiner BL, Phillips SF. Increased eosinophil granule proteins in gut lavage fluid from patients with inflammatory bowel disease. *Mayo Clin Proc*. Feb 1997;72(2):117-123.
71. Bischoff SC, Grabowsky J, Manns MP. Quantification of inflammatory mediators in stool samples of patients with inflammatory bowel disorders and controls. *Dig Dis Sci*. Feb 1997;42(2):394-403.
72. Saitoh O, Kojima K, Sugi K, et al. Fecal eosinophil granule-derived proteins reflect disease activity in inflammatory bowel disease. *Am J Gastroenterol*. Dec 1999;94(12):3513-3520.
73. Liu LX, Chi J, Upton MP, Ash LR. Eosinophilic colitis associated with larvae of the pinworm *Enterobius vermicularis*. *Lancet*. Aug 12 1995;346(8972):410-412.
74. Bischoff SC, Mayer J, Nguyen QT, Stolte M, Manns MP. Immunohistological assessment of intestinal eosinophil activation in patients with eosinophilic gastroenteritis and inflammatory bowel disease. *Am J Gastroenterol*. 12/1999 1999;94(12):3521-3529.
75. Clouse RE, Alpers DH, Hockenbery DM, DeSchryver-Kecskemeti K. Pericrypt eosinophilic enterocolitis and chronic diarrhea. *Gastroenterology*. Jul 1992;103(1):168-176.
76. Bemark M, Boysen P, Lycke NY. Induction of gut IgA production through T cell-dependent and T cell-independent pathways. *Annals of the New York Academy of Sciences*. Jan 2012;1247:97-116.

77. Faria AM, Gomes-Santos AC, Goncalves JL, et al. Food components and the immune system: from tonic agents to allergens. *Frontiers in immunology*. 2013;4:102.
78. Corthesy B. Secretory immunoglobulin A: well beyond immune exclusion at mucosal surfaces. *Immunopharmacology and immunotoxicology*. Jun 2009;31(2):174-179.
79. Everett ML, Palestrant D, Miller SE, Bollinger RR, Parker W. Immune exclusion and immune inclusion: A new model of host-bacterial interactions in the gut. *Clinical and Applied Immunology Reviews* March 15, 2004;4:321-332.
80. Bollinger RR EM, Palestrant D. Human secretory immunoglobulin A may contribute to biofilm formation in the gut. *Immunology*. 2003;109(4):580-587.
81. Shroff Ke MKCJ. Commensal enteric bacteria engender a self-limiting humoral mucosal immune response while permanently colonizing the gut. *Infect Immun*. . 1995;63:3904-3913.
82. Neutra MR, Kraehenbuhl JP. Mucosal immunization via M cells for production of protective secretory IgA antibodies. *The American journal of tropical medicine and hygiene*. 1994;50(5 Suppl):10-13.
83. Bollinger RR, Everett ML, Palestrant D, Love SD, Lin SS, Parker W. Human secretory immunoglobulin A may contribute to biofilm formation in the gut. *Immunology*. Aug 2003;109(4):580-587.
84. Martins FS, Rodrigues AC, Tiago FC, et al. Saccharomyces cerevisiae strain 905 reduces the translocation of Salmonella enterica serotype Typhimurium and stimulates the immune system in gnotobiotic and conventional mice. *Journal of medical microbiology*. Mar 2007;56(Pt 3):352-359.
85. Frossard CP, Hauser C, Eigenmann PA. Antigen-specific secretory IgA antibodies in the gut are decreased in a mouse model of food allergy. *The Journal of allergy and clinical immunology*. Aug 2004;114(2):377-382.
86. Maruyama K, Hida M, Kohgo T, Fukunaga Y. Changes in salivary and fecal secretory IgA in infants under different feeding regimens. *Pediatr Int*. Jun 2009;51(3):342-345.
87. Hauk PJ. The role of food allergy in atopic dermatitis. *Curr Allergy Asthma Rep*. May 2008;8(3):188-194.
88. Jarvinen KM, Konstantinou GN, Pilpil M, et al. Intestinal permeability in children with food allergy on specific elimination diets. *Pediatr Allergy Immunol*. Sep 2013;24(6):589-595.
89. Perrier C, Corthesy B. Gut permeability and food allergies. *Clinical and experimental allergy : journal of the British Society for Allergy and Clinical Immunology*. Jan 2011;41(1):20-28.
90. Yu LC. The epithelial gatekeeper against food allergy. *Pediatrics and neonatology*. Dec 2009;50(6):247-254.
91. Scalabrin DM, Mitmesser SH, Welling GW, et al. New prebiotic blend of polydextrose and galacto-oligosaccharides has a bifidogenic effect in young infants. *Journal of pediatric gastroenterology and nutrition*. Mar 2012;54(3):343-352.
92. Mohan R, Koebnick C, Schildt J, Mueller M, Radke M, Blaut M. Effects of Bifidobacterium lactis Bb12 supplementation on body weight, fecal pH, acetate, lactate, calprotectin, and IgA in preterm infants. *Pediatric research*. Oct 2008;64(4):418-422.
93. Bakker-Zierikzee AM, Tol EA, Kroes H, Alles MS, Kok FJ, Bindels JG. Faecal SIgA secretion in infants fed on pre- or probiotic infant formula. *Pediatr Allergy Immunol*. Mar 2006;17(2):134-140.
94. Saarinen KM, Sarnesto A, Savilahti E. Markers of inflammation in the feces of infants with cow's milk allergy. *Pediatric allergy and immunology : official publication of the European Society of Pediatric Allergy and Immunology*. Jun 2002;13(3):188-194.
95. Glocker E, Ehl S, Grimbacher B. Common variable immunodeficiency in children. *Current opinion in pediatrics*. Dec 2007;19(6):685-692.
96. Fagarasan S, Honjo T. Intestinal IgA synthesis: regulation of front-line body defences. *Nat Rev Immunol*. Jan 2003;3(1):63-72.
97. Martins CA, Fonteles MG, Barrett LJ, Guerrant RL. Correlation of lactoferrin with neutrophilic inflammation in body fluids. *Clin Diagn Lab Immunol*. 11/1995 1995;2(6):763-765.
98. Conneely OM. Antiinflammatory activities of lactoferrin. *Journal of the American College of Nutrition*. 10/2001 2001;20(5 Suppl):389S-395S.
99. Fine KD, Ogunji F, George J, Niehaus MD, Guerrant RL. Utility of a rapid fecal latex agglutination test detecting the neutrophil protein, lactoferrin, for diagnosing inflammatory causes of chronic diarrhea. *The American journal of gastroenterology*. Aug 1998;93(8):1300-1305.
100. Sidhu R, Sanders DS, Wilson P, Foye L, Morley S, McAlindon ME. Faecal lactoferrin, capsule endoscopy and Crohn's disease. Is there a three way relationship? A pilot study. *J Gastrointest Liver Dis*. Sep 2010;19(3):257-260.
101. van Langenberg DR, Gearry RB, Wong HL, Ward M, Gibson PR. The potential value of faecal lactoferrin as a screening test in hospitalized patients with diarrhoea. *Internal medicine journal*. Dec 2010;40(12):819-827.
102. Sidhu R, Wilson P, Wright A, et al. Faecal lactoferrin--a novel test to differentiate between the irritable and inflamed bowel? *Alimentary pharmacology & therapeutics*. Jun 2010;31(12):1365-1370.
103. Vieira A, Fang CB, Rolim EG, et al. Inflammatory bowel disease activity assessed by fecal calprotectin and lactoferrin: correlation with laboratory parameters, clinical, endoscopic and histological indexes. *BMC.Res.Notes*. 2009 2009;2:221.
104. Kim HS, Lim JH, Park H, Lee SI. Increased immunoendocrine cells in intestinal mucosa of postinfectious irritable bowel syndrome patients 3 years after acute Shigella infection--an observation in a small case control study. *Yonsei Med J*. Jan 2010;51(1):45-51.
105. Rupa P, Mine Y. Recent Advances in the Role of Probiotics in Human Inflammation and Gut Health. *Journal of agricultural and food chemistry*. Aug 16 2012.
106. Quigley EM. Therapies aimed at the gut microbiota and inflammation: antibiotics, prebiotics, probiotics, synbiotics, anti-inflammatory therapies. *Gastroenterol Clin North Am*. Mar 2011;40(1):207-222.
107. Fukushima Y, Kawata Y, Hara H, Terada A, Mitsuoka T. Effect of a probiotic formula on intestinal immunoglobulin A production in healthy children. *International journal of food microbiology*. Jun 30 1998;42(1-2):39-44.
108. Khosraviani K, Weir HP, Hamilton P, Moorehead J, Williamson K. Effect of folate supplementation on mucosal cell proliferation in high risk patients for colon cancer. *Gut*. Aug 2002;51(2):195-199.
109. Vermeulen MA, de Jong J, Vaessen MJ, van Leeuwen PA, Houdijk AP. Glutamate reduces experimental intestinal hyperpermeability and facilitates glutamine support of gut integrity. *World J Gastroenterol*. Mar 28 2011;17(12):1569-1573.
110. Coeffier M, Claeysens S, Hecketsweiler B, Lavoigne A, Ducrotte P, Dechelotte P. Enteral glutamine stimulates protein synthesis and decreases ubiquitin mRNA level in human gut mucosa. *American journal of physiology. Gastrointestinal and liver physiology*. Aug 2003;285(2):G266-273.
111. Lewis K LF, Phan V, Söderholm JD, Sherman PM, McKay DM. Enhanced translocation of bacteria across metabolically stressed epithelia is reduced by butyrate. *Inflamm Bowel Dis*. . 2010;16(7):1138-1148.
112. Abrahao V. Nourishing the dysfunctional gut and whey protein. *Current opinion in clinical nutrition and metabolic care*. Sep 2012;15(5):480-484.
113. Hawrelak JA MS. Effects of two natural medicine formulations on irritable bowel syndrome symptoms: a pilot study. *J Altern Complement Med*. 2010;Oct;16(10):1065-1071.



114. Braga M. Perioperative immunonutrition and gut function. *Current opinion in clinical nutrition and metabolic care*. Sep 2012;15(5):485-488.
115. Taylor RA, Leonard MC. Curcumin for inflammatory bowel disease: a review of human studies. *Alternative medicine review : a journal of clinical therapeutic*. Jun 2011;16(2):152-156.
116. Carroccio A, Mansueto P, Iacono G, et al. Non-celiac wheat sensitivity diagnosed by double-blind placebo-controlled challenge: exploring a new clinical entity. *Am J Gastroenterol*. Dec 2012;107(12):1898-1906; quiz 1907.
117. Guo H, Jiang T, Wang J, Chang Y, Guo H, Zhang W. The value of eliminating foods according to food-specific immunoglobulin G antibodies in irritable bowel syndrome with diarrhoea. *The Journal of international medical research*. 2012;40(1):204-210.
118. Zar S, Mincher L, Benson MJ, Kumar D. Food-specific IgG4 antibody-guided exclusion diet improves symptoms and rectal compliance in irritable bowel syndrome. *Scand J Gastroenterol*. Jul 2005;40(7):800-807.
119. Lason A, Simren M, Stotzer PO, Isaksson S, Ohman L, Strid H. Fecal calprotectin levels predict the clinical course in patients with new onset of ulcerative colitis. *Inflammatory bowel diseases*. Mar 2013;19(3):576-581.
120. Thomas LV, Ockhuizen T. New insights into the impact of the intestinal microbiota on health and disease: a symposium report. *The British journal of nutrition*. Jan 2012;107 Suppl 1:S1-13.
121. Young VB. The intestinal microbiota in health and disease. *Current opinion in gastroenterology*. Jan 2012;28(1):63-69.
122. Clarke G, Stilling RM, Kennedy PJ, Stanton C, Cryan JF, Dinan TG. Gut Microbiota: The Neglected Endocrine Organ. *Molecular endocrinology (Baltimore, Md.)*. Jun 3 2014;me20141108.
123. Dulal S, Keku TO. Gut microbiome and colorectal adenomas. *Cancer journal (Sudbury, Mass.)*. May-Jun 2014;20(3):225-231.
124. Ebel B, Lemetais G, Beney L, et al. Impact of probiotics on risk factors for cardiovascular diseases. A review. *Crit Rev Food Sci Nutr*. 2014;54(2):175-189.
125. Ridaura VK, Faith JJ, Rey FE, et al. Gut microbiota from twins discordant for obesity modulate metabolism in mice. *Science (New York, N.Y.)*. Sep 6 2013;341(6150):1241-1244.
126. Xiao S, Zhao L. Gut microbiota-based translational biomarkers to prevent metabolic syndrome via nutritional modulation. *FEMS Microbiol Ecol*. Feb 2014;87(2):303-314.
127. Wong JM, de Souza R, Kendall CW, Emam A, Jenkins DJ. Colonic health: fermentation and short chain fatty acids. *Journal of clinical gastroenterology*. Mar 2006;40(3):235-243.
128. Mortensen PB, Clausen MR. Short-chain fatty acids in the human colon: relation to gastrointestinal health and disease. *Scandinavian journal of gastroenterology. Supplement*. 1996;216:132-148.
129. Velazquez OC, Lederer HM, Rombeau JL. Butyrate and the colonocyte. Production, absorption, metabolism, and therapeutic implications. *Advances in experimental medicine and biology*. 1997;427:123-134.
130. Cummings JH, Pomare EW, Branch WJ, Naylor CP, Macfarlane GT. Short chain fatty acids in human large intestine, portal, hepatic and venous blood. *Gut*. Oct 1987;28(10):1221-1227.
131. Tibble J, Teahon K, Thjodleifsson B, et al. A simple method for assessing intestinal inflammation in Crohn's disease. *Gut*. Oct 2000;47(4):506-513.
132. Bunn SK, Bisset WM, Main MJ, Golden BE. Fecal calprotectin as a measure of disease activity in childhood inflammatory bowel disease. *J Pediatr Gastroenterol Nutr*. Feb 2001;32(2):171-177.
133. Walker AW, Lawley TD. Therapeutic modulation of intestinal dysbiosis. *Pharmacological research : the official journal of the Italian Pharmacological Society*. Sep 24 2012.
134. Daulatzai MA. Chronic functional bowel syndrome enhances gut-brain axis dysfunction, neuroinflammation, cognitive impairment, and vulnerability to dementia. *Neurochemical research*. Apr 2014;39(4):624-644.
135. Goepf J, Fowler E, McBride T, Landis D. Frequency of abnormal fecal biomarkers in irritable bowel syndrome. *Global advances in health and medicine : improving healthcare outcomes worldwide*. May 2014;3(3):9-15.
136. Lopez-Siles M, Martinez-Medina M, Busquets D, et al. Mucosa-associated Faecalibacterium prausnitzii and Escherichia coli co-abundance can distinguish Irritable Bowel Syndrome and Inflammatory Bowel Disease phenotypes. *International journal of medical microbiology : IJMM*. May 2014;304(3-4):464-475.
137. Klurfeld DM. Fiber and cancer protection--mechanisms. *Adv Exp Med Biol*. 1997;427:249-257.
138. MacLennan R, Jensen OM, Mosbech J, Vuori H. Diet, transit time, stool weight, and colon cancer in two Scandinavian populations. *Am J Clin Nutr*. Oct 1978;31(10 Suppl):S239-S242.
139. MacLennan R, Jensen OM. Dietary fibre, transit-time, faecal bacteria, steroids, and colon cancer in two Scandinavian populations. Report from the International Agency for Research on Cancer Intestinal Microecology Group. *Lancet*. Jul 30 1977;2(8031):207-211.
140. Walker AR. Colon cancer and diet, with special reference to intakes of fat and fiber. *Am J Clin Nutr*. Dec 1976;29(12):1417-1426.
141. Modan B, Barel V, Lubin F, Modan M. Dietary factors and cancer in Israel. *Cancer research*. Nov 1975;35(11 Pt. 2):3503-3506.
142. Howe GR, Benito E, Castelleto R, et al. Dietary intake of fiber and decreased risk of cancers of the colon and rectum: evidence from the combined analysis of 13 case-control studies. *J Natl Cancer Inst*. Dec 16 1992;84(24):1887-1896.
143. Reddy BS, Ekelund G, Bohe M, Engle A, Domellof L. Metabolic epidemiology of colon cancer: dietary pattern and fecal sterol concentrations of three populations. *Nutrition and cancer*. 1983;5(1):34-40.
144. Tollefson L. The use of epidemiology, scientific data, and regulatory authority to determine risk factors in cancer of some organs of the digestive system. 4. Colon cancer. *Regul Toxicol Pharmacol*. Mar 1986;6(1):24-54.
145. Doheny K. High-Fiber Diet Linked to Lower Colon Cancer Risk: Research Shows Fiber From Whole Grains and Cereals Associated With Reduced Colon Cancer Risk. Nov. 10, 2011 <http://www.webmd.com/colorectal-cancer/news/20111110/high-fiber-diet-linked-to-lower-colon-cancer-risk>.
146. McOrist AL, Miller RB, Bird AR, et al. Fecal butyrate levels vary widely among individuals but are usually increased by a diet high in resistant starch. *The Journal of nutrition*. May 2011;141(5):883-889.
147. Bjarnason I, Sherwood R. Fecal calprotectin: a significant step in the noninvasive assessment of intestinal inflammation. *J Pediatr Gastroenterol Nutr*. Jul 2001;33(1):11-13.
148. Reddy BS, Hedges AR, Laakso K, Wynder EL. Metabolic epidemiology of large bowel cancer: fecal bulk and constituents of high-risk North American and low-risk Finnish population. *Cancer*. Dec 1978;42(6):2832-2838.
149. Britton RA, Young VB. Role of the Intestinal Microbiota in Resistance to Colonization by Clostridium difficile. *Gastroenterology*. Feb 3 2014.
150. Buffie CG, Pamer EG. Microbiota-mediated colonization resistance against intestinal pathogens. *Nature reviews*. Nov 2013;13(11):790-801.
151. Brown K, DeCoffe D, Molcan E, Gibson DL. Diet-induced dysbiosis of the intestinal microbiota and the effects on immunity and disease. *Nutrients*. Aug 2012;4(8):1095-1119.

152. Chassard C, Dapoigny M, Scott KP, et al. Functional dysbiosis within the gut microbiota of patients with constipated-irritable bowel syndrome. *Alimentary pharmacology & therapeutics*. Apr 2012;35(7):828-838.
153. Comito D, Romano C. Dysbiosis in the pathogenesis of pediatric inflammatory bowel diseases. *Int J Inflam*. 2012;2012:687143.
154. Craven M, Egan CE, Dowd SE, et al. Inflammation drives dysbiosis and bacterial invasion in murine models of ileal Crohn's disease. *PLoS One*. 2012;7(7):e41594.
155. Duboc H, Rajca S, Rainteau D, et al. Connecting dysbiosis, bile-acid dysmetabolism and gut inflammation in inflammatory bowel diseases. *Gut*. Sep 19 2012.
156. Fujiyama Y, Andoh A. [Dysbiosis in inflammatory bowel disease]. *Nihon Rinsho*. Feb 2012;70 Suppl 1:79-84.
157. Joossens M, et al. . Dysbiosis of the faecal microbiota in patients with Crohn's disease and thier unaffected relatives. *Gut*. 2012;60:631-637.
158. Kang S, Denman SE, Morrison M, et al. Dysbiosis of fecal microbiota in Crohn's disease patients as revealed by a custom phylogenetic microarray. *Inflamm Bowel Dis*. Dec 2010;16(12):2034-2042.
159. Kaur N, Chen CC, Luther J, Kao JY. Intestinal dysbiosis in inflammatory bowel disease. *Gut Microbes*. Jul-Aug 2011;2(4):211-216.
160. Ogilvie LA, Jones BV. Dysbiosis modulates capacity for bile acid modification in the gut microbiomes of patients with inflammatory bowel disease: a mechanism and marker of disease? *Gut*. Nov 2012;61(11):1642-1643.
161. Packey CD, Sartor RB. Commensal bacteria, traditional and opportunistic pathogens, dysbiosis and bacterial killing in inflammatory bowel diseases. *Curr Opin Infect Dis*. Jun 2009;22(3):292-301.
162. Sun L, Nava GM, Stappenbeck TS. Host genetic susceptibility, dysbiosis, and viral triggers in inflammatory bowel disease. *Current opinion in gastroenterology*. Jul 2011;27(4):321-327.
163. Walker AW, Lawley TD. Therapeutic modulation of intestinal dysbiosis. *Pharmacol Res*. Sep 24 2012.
164. Duncan SH, Belenguer A, Holtrop G, Johnstone AM, Flint HJ, Loble GE. Reduced dietary intake of carbohydrates by obese subjects results in decreased concentrations of butyrate and butyrate-producing bacteria in feces. *Applied and environmental microbiology*. Feb 2007;73(4):1073-1078.
165. Zimmer J, Lange B, Frick JS, et al. A vegan or vegetarian diet substantially alters the human colonic faecal microbiota. *European journal of clinical nutrition*. Jan 2012;66(1):53-60.
166. Sanz Y. Effects of a gluten-free diet on gut microbiota and immune function in healthy adult humans. *Gut Microbes*. May-Jun 2010;1(3):135-137.
167. De Palma G, Nadal I, Collado MC, Sanz Y. Effects of a gluten-free diet on gut microbiota and immune function in healthy adult human subjects. *The British journal of nutrition*. Oct 2009;102(8):1154-1160.
168. Nadal I, Donat E, Ribes-Koninckx C, Calabuig M, Sanz Y. Imbalance in the composition of the duodenal microbiota of children with coeliac disease. *Journal of medical microbiology*. Dec 2007;56(Pt 12):1669-1674.
169. Collado MC, Donat E, Ribes-Koninckx C, Calabuig M, Sanz Y. Specific duodenal and faecal bacterial groups associated with paediatric coeliac disease. *Journal of clinical pathology*. Mar 2009;62(3):264-269.
170. Collado MC, Donat E, Ribes-Koninckx C, Calabuig M, Sanz Y. Imbalances in faecal and duodenal Bifidobacterium species composition in active and non-active coeliac disease. *BMC Microbiol*. 2008;8:232.
171. David LA, Maurice CF, Carmody RN, et al. Diet rapidly and reproducibly alters the human gut microbiome. *Nature*. Jan 23 2014;505(7484):559-563.
172. Le Chatelier E, Nielsen T, Qin J, et al. Richness of human gut microbiome correlates with metabolic markers. *Nature*. Aug 29 2013;500(7464):541-546.
173. Noor SO, Ridgway K, Scovell L, et al. Ulcerative colitis and irritable bowel patients exhibit distinct abnormalities of the gut microbiota. *BMC gastroenterology*. 2010;10:134.
174. Benjamin JL, Hedin CR, Koutsoumpas A, et al. Smokers with active Crohn's disease have a clinically relevant dysbiosis of the gastrointestinal microbiota. *Inflammatory bowel diseases*. Jun 2012;18(6):1092-1100.
175. Scanlan PD, Shanahan F, O'Mahony C, Marchesi JR. Culture-independent analyses of temporal variation of the dominant fecal microbiota and targeted bacterial subgroups in Crohn's disease. *Journal of clinical microbiology*. Nov 2006;44(11):3980-3988.
176. Andoh A, Sakata S, Koizumi Y, Mitsuyama K, Fujiyama Y, Benno Y. Terminal restriction fragment length polymorphism analysis of the diversity of fecal microbiota in patients with ulcerative colitis. *Inflammatory bowel diseases*. Aug 2007;13(8):955-962.
177. Swidsinski A, Ladhoff A, Pernthaler A, et al. Mucosal flora in inflammatory bowel disease. *Gastroenterology*. Jan 2002;122(1):44-54.
178. Vahtovuo J, Munukka E, Korkeamaki M, Luukkainen R, Toivanen P. Fecal microbiota in early rheumatoid arthritis. *The Journal of rheumatology*. Aug 2008;35(8):1500-1505.
179. Collado MC, Isolauri E, Laitinen K, Salminen S. Distinct composition of gut microbiota during pregnancy in overweight and normal-weight women. *The American journal of clinical nutrition*. Oct 2008;88(4):894-899.
180. Schwiertz A, Taras D, Schafer K, et al. Microbiota and SCFA in lean and overweight healthy subjects. *Obesity (Silver Spring)*. Jan 2010;18(1):190-195.
181. Patil DP, Dhotre DP, Chavan SG, et al. Molecular analysis of gut microbiota in obesity among Indian individuals. *J Biosci*. Sep 2012;37(4):647-657.
182. Furet JB, Kong LC, Tap J, et al. Differential adaptation of human gut microbiota to bariatric surgery-induced weight loss: links with metabolic and low-grade inflammation markers. *Diabetes*. Dec 2010;59(12):3049-3057.
183. Larsen N, Vogensen FK, van den Berg FW, et al. Gut microbiota in human adults with type 2 diabetes differs from non-diabetic adults. *PLoS One*. 2010;5(2):e9085.
184. Zhang X, Shen D, Fang Z, et al. Human gut microbiota changes reveal the progression of glucose intolerance. *PLoS One*. 2013;8(8):e71108.
185. Samb-Ba B, Mazonot C, Gassama-Sow A, et al. MALDI-TOF Identification of the Human Gut Microbiome in People with and without Diarrhea in Senegal. *PloS one*. 2014;9(5):e87419.
186. Bervoets L, Van Hoorenbeeck K, Kortleven I, et al. Differences in gut microbiota composition between obese and lean children: a cross-sectional study. *Gut Pathog*. 2013;5(1):10.
187. Wu X, Ma C, Han L, et al. Molecular characterisation of the faecal microbiota in patients with type II diabetes. *Current microbiology*. Jul 2010;61(1):69-78.
188. Finegold SM, Dowd SE, Gontcharova V, et al. Pyrosequencing study of fecal microflora of autistic and control children. *Anaerobe*. Aug 2010;16(4):444-453.
189. Sanchez E, De Palma G, Capilla A, et al. Influence of environmental and genetic factors linked to celiac disease risk on infant gut colonization by Bacteroides species. *Applied and environmental microbiology*. Aug 2011;77(15):5316-5323.
190. Sanchez E, Donat E, Ribes-Koninckx C, Calabuig M, Sanz Y. Intestinal Bacteroides species associated with coeliac disease. *Journal of clinical pathology*. Dec 2010;63(12):1105-1111.



191. Wylie KM, Truty RM, Sharpton TJ, et al. Novel bacterial taxa in the human microbiome. *PLoS One*. 2012;7(6):e35294.
192. Ubeda C, Bucci V, Caballero S, et al. Intestinal microbiota containing *Barnesiella* species cures vancomycin-resistant *Enterococcus faecium* colonization. *Infect Immun*. Mar 2013;81(3):965-973.
193. Morgan XC, Tickle TL, Sokol H, et al. Dysfunction of the intestinal microbiome in inflammatory bowel disease and treatment. *Genome biology*. 2012;13(9):R79.
194. Swidsinski A, Weber J, Loening-Baucke V, Hale LP, Lochs H. Spatial organization and composition of the mucosal flora in patients with inflammatory bowel disease. *Journal of clinical microbiology*. Jul 2005;43(7):3380-3389.
195. Scher JU, Szczesnak A, Longman RS, et al. Expansion of intestinal *Prevotella copri* correlates with enhanced susceptibility to arthritis. *eLife*. 2013;2:e01202.
196. Kang DW, Park JG, Ilhan ZE, et al. Reduced incidence of *Prevotella* and other fermenters in intestinal microflora of autistic children. *PLoS One*. 2013;8(7):e68322.
197. Claesson MJ, Jeffery IB, Conde S, et al. Gut microbiota composition correlates with diet and health in the elderly. *Nature*. Aug 9 2012;488(7410):178-184.
198. Zhu L, Baker SS, Gill C, et al. Characterization of gut microbiomes in nonalcoholic steatohepatitis (NASH) patients: a connection between endogenous alcohol and NASH. *Hepatology (Baltimore, Md.)*. Feb 2013;57(2):601-609.
199. Zupancic ML, Cantarel BL, Liu Z, et al. Analysis of the gut microbiota in the old order Amish and its relation to the metabolic syndrome. *PLoS One*. 2012;7(8):e43052.
200. Rajilic-Stojanovic M, Biagi E, Heilig HG, et al. Global and deep molecular analysis of microbiota signatures in fecal samples from patients with irritable bowel syndrome. *Gastroenterology*. Nov 2011;141(5):1792-1801.
201. Song Y, Liu C, Finegold SM. Real-time PCR quantitation of clostridia in feces of autistic children. *Applied and environmental microbiology*. Nov 2004;70(11):6459-6465.
202. Finegold SM, Molitoris D, Song Y, et al. Gastrointestinal microflora studies in late-onset autism. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America*. Sep 1 2002;35(Suppl 1):S6-S16.
203. Parracho HM, Bingham MO, Gibson GR, McCartney AL. Differences between the gut microflora of children with autistic spectrum disorders and that of healthy children. *Journal of medical microbiology*. Oct 2005;54(Pt 10):987-991.
204. Sandler RH, Finegold SM, Bolte ER, et al. Short-term benefit from oral vancomycin treatment of regressive-onset autism. *Journal of child neurology*. Jul 2000;15(7):429-435.
205. Fava F, Danese S. Intestinal microbiota in inflammatory bowel disease: friend of foe? *World J Gastroenterol*. Feb 7 2011;17(5):557-566.
206. Williams BL, Hornig M, Buie T, et al. Impaired carbohydrate digestion and transport and mucosal dysbiosis in the intestines of children with autism and gastrointestinal disturbances. *PLoS One*. 2011;6(9):e24585.
207. Malinen E, Krogius-Kurikka L, Lyra A, et al. Association of symptoms with gastrointestinal microbiota in irritable bowel syndrome. *World J Gastroenterol*. Sep 28 2010;16(36):4532-4540.
208. Miquel S, Martin R, Rossi O, et al. *Faecalibacterium prausnitzii* and human intestinal health. *Current opinion in microbiology*. Jun 2013;16(3):255-261.
209. Sokol H, Seksik P, Furet JP, et al. Low counts of *Faecalibacterium prausnitzii* in colitis microbiota. *Inflammatory bowel diseases*. Aug 2009;15(8):1183-1189.
210. Berry D, Reinisch W. Intestinal microbiota: a source of novel biomarkers in inflammatory bowel diseases? *Best practice & research. Clinical gastroenterology*. Feb 2013;27(1):47-58.
211. Sokol H, Pigneur B, Watterlot L, et al. *Faecalibacterium prausnitzii* is an anti-inflammatory commensal bacterium identified by gut microbiota analysis of Crohn disease patients. *Proceedings of the National Academy of Sciences of the United States of America*. Oct 28 2008;105(43):16731-16736.
212. Joossens M, Huys G, Cnockaert M, et al. Dysbiosis of the faecal microbiota in patients with Crohn's disease and their unaffected relatives. *Gut*. May 2011;60(5):631-637.
213. Schwartz A, Jacobi M, Frick JS, Richter M, Rusch K, Kohler H. Microbiota in pediatric inflammatory bowel disease. *The Journal of pediatrics*. Aug 2010;157(2):240-244 e241.
214. Frank DN, Robertson CE, Hamm CM, et al. Disease phenotype and genotype are associated with shifts in intestinal-associated microbiota in inflammatory bowel diseases. *Inflammatory bowel diseases*. Jan 2011;17(1):179-184.
215. Swidsinski A, Loening-Baucke V, Vaneechoutte M, Doerffel Y. Active Crohn's disease and ulcerative colitis can be specifically diagnosed and monitored based on the biostructure of the fecal flora. *Inflammatory bowel diseases*. Feb 2008;14(2):147-161.
216. Machiels K, Joossens M, Sabino J, et al. A decrease of the butyrate-producing species *Roseburia hominis* and *Faecalibacterium prausnitzii* defines dysbiosis in patients with ulcerative colitis. *Gut*. Sep 10 2013.
217. Maccaferri S, Candela M, Turroni S, et al. IBS-associated phylogenetic unbalances of the intestinal microbiota are not reverted by probiotic supplementation. *Gut microbes*. Sep-Oct 2012;3(5):406-413.
218. Tana C, Umesaki Y, Imaoka A, Handa T, Kanazawa M, Fukudo S. Altered profiles of intestinal microbiota and organic acids may be the origin of symptoms in irritable bowel syndrome. *Neurogastroenterology and motility : the official journal of the European Gastrointestinal Motility Society*. May 2010;22(5):512-519, e114-515.
219. Rigsbee L, Agans R, Shankar V, et al. Quantitative profiling of gut microbiota of children with diarrhea-predominant irritable bowel syndrome. *The American journal of gastroenterology*. Nov 2012;107(11):1740-1751.
220. Balsari A, Ceccarelli A, Dubini F, Fesce E, Poli G. The fecal microbial population in the irritable bowel syndrome. *Microbiologica*. Jul 1982;5(3):185-194.
221. Malinen E, Rinttila T, Kajander K, et al. Analysis of the fecal microbiota of irritable bowel syndrome patients and healthy controls with real-time PCR. *The American journal of gastroenterology*. Feb 2005;100(2):373-382.
222. Krogius-Kurikka L, Lyra A, Malinen E, et al. Microbial community analysis reveals high level phylogenetic alterations in the overall gastrointestinal microbiota of diarrhoea-predominant irritable bowel syndrome sufferers. *BMC gastroenterology*. 2009;9:95.
223. Armougom F, Henry M, Vialettes B, Raccach D, Raoult D. Monitoring bacterial community of human gut microbiota reveals an increase in *Lactobacillus* in obese patients and *Methanogens* in anorexic patients. *PLoS One*. 2009;4(9):e7125.
224. Chassard C, Dapoigny M, Scott KP, et al. Functional dysbiosis within the gut microbiota of patients with constipated-irritable bowel syndrome. *Alimentary pharmacology & therapeutics*. Apr 2012;35(7):828-838.
225. Png CW, Linden SK, Gilshenan KS, et al. Mucolytic bacteria with increased prevalence in IBD mucosa augment in vitro utilization of mucin by other bacteria. *The American journal of gastroenterology*. Nov 2010;105(11):2420-2428.
226. Joossens M, Huys G, Cnockaert M, et al. Dysbiosis of the faecal microbiota in patients with Crohn's disease and their unaffected relatives. *Gut*. 2011;60(5):631-637.

227. Willing BP, Dicksved J, Halfvarson J, et al. A pyrosequencing study in twins shows that gastrointestinal microbial profiles vary with inflammatory bowel disease phenotypes. *Gastroenterology*. Dec 2010;139(6):1844-1854 e1841.
228. Lyra A, Rinttila T, Nikkila J, et al. Diarrhoea-predominant irritable bowel syndrome distinguishable by 16S rRNA gene phylotype quantification. *World J Gastroenterol*. Dec 21 2009;15(47):5936-5945.
229. Rossi M, Amaretti A, Raimondi S. Folate production by probiotic bacteria. *Nutrients*. Jan 2011;3(1):118-134.
230. Kerckhoffs AP, Samsom M, van der Rest ME, et al. Lower Bifidobacteria counts in both duodenal mucosa-associated and fecal microbiota in irritable bowel syndrome patients. *World J Gastroenterol*. Jun 21 2009;15(23):2887-2892.
231. Si JM, Yu YC, Fan YJ, Chen SJ. Intestinal microecology and quality of life in irritable bowel syndrome patients. *World J Gastroenterol*. Jun 15 2004;10(12):1802-1805.
232. Xu X, Hui H, Cai D. [Differences in fecal Bifidobacterium species between patients with type 2 diabetes and healthy individuals]. *Nan fang yi ke da xue xue bao = Journal of Southern Medical University*. Apr 2012;32(4):531-533, 564.
233. Le KA, Li Y, Xu X, et al. Alterations in fecal Lactobacillus and Bifidobacterium species in type 2 diabetic patients in Southern China population. *Frontiers in physiology*. 2012;3:496.
234. Sekirov I, Russell SL, Antunes LC, Finlay BB. Gut microbiota in health and disease. *Physiological reviews*. Jul 2010;90(3):859-904.
235. Wang L, Christophersen CT, Sorich MJ, Gerber JP, Anglely MT, Conlon MA. Low relative abundances of the mucolytic bacterium Akkermansia muciniphila and Bifidobacterium spp. in feces of children with autism. *Applied and environmental microbiology*. Sep 2011;77(18):6718-6721.
236. Turnbaugh PJ, Hamady M, Yatsunenko T, et al. A core gut microbiome in obese and lean twins. *Nature*. Jan 22 2009;457(7228):480-484.
237. Kalliomaki M, Collado MC, Salminen S, Isolauri E. Early differences in fecal microbiota composition in children may predict overweight. *The American journal of clinical nutrition*. Mar 2008;87(3):534-538.
238. Santacruz A, Marcos A, Warnberg J, et al. Interplay between weight loss and gut microbiota composition in overweight adolescents. *Obesity (Silver Spring)*. Oct 2009;17(10):1906-1915.
239. Kassinen A, Krogius-Kurikka L, Makivuokko H, et al. The fecal microbiota of irritable bowel syndrome patients differs significantly from that of healthy subjects. *Gastroenterology*. Jul 2007;133(1):24-33.
240. Loubinoux J, Bronowicki JP, Pereira IA, Mouguel JL, Faou AE. Sulfate-reducing bacteria in human feces and their association with inflammatory bowel diseases. *FEMS microbiology ecology*. May 1 2002;40(2):107-112.
241. Roediger WE, Moore J, Babbage W. Colonic sulfide in pathogenesis and treatment of ulcerative colitis. *Digestive diseases and sciences*. Aug 1997;42(8):1571-1579.
242. Finegold SM, Downes J, Summanen PH. Microbiology of regressive autism. *Anaerobe*. Apr 2012;18(2):260-262.
243. Karlsson CL, Onnerfalt J, Xu J, Molin G, Ahrne S, Thorngren-Jerneck K. The microbiota of the gut in preschool children with normal and excessive body weight. *Obesity (Silver Spring)*. Nov 2012;20(11):2257-2261.
244. Tiihonen K, Ouwehand AC, Rautonen N. Effect of overweight on gastrointestinal microbiology and immunology: correlation with blood biomarkers. *Br J Nutr*. Apr 2010;103(7):1070-1078.
245. Duboc H, Rajca S, Rainteau D, et al. Connecting dysbiosis, bile-acid dysmetabolism and gut inflammation in inflammatory bowel diseases. *Gut*. Apr 2013;62(4):531-539.
246. Iebba V, Aloï M, Civitelli F, Cucchiara S. Gut microbiota and pediatric disease. *Digestive diseases (Basel, Switzerland)*. 2011;29(6):531-539.
247. Comito D, Romano C. Dysbiosis in the pathogenesis of pediatric inflammatory bowel diseases. *International journal of inflammation*. 2012;2012:687143.
248. Duboc H, Rainteau D, Rajca S, et al. Increase in fecal primary bile acids and dysbiosis in patients with diarrhea-predominant irritable bowel syndrome. *Neurogastroenterology and motility : the official journal of the European Gastrointestinal Motility Society*. Jun 2012;24(6):513-520, e246-517.
249. Santacruz A, Collado MC, Garcia-Valdes L, et al. Gut microbiota composition is associated with body weight, weight gain and biochemical parameters in pregnant women. *Br J Nutr*. Jul 2010;104(1):83-92.
250. Kaufman DW, Kelly JP, Curhan GC, et al. Oxalobacter formigenes may reduce the risk of calcium oxalate kidney stones. *Journal of the American Society of Nephrology : JASN*. Jun 2008;19(6):1197-1203.
251. Knight J, Deora R, Assimos DG, Holmes RP. The genetic composition of Oxalobacter formigenes and its relationship to colonization and calcium oxalate stone disease. *Urolithiasis*. Jun 2013;41(3):187-196.
252. Siener R, Bangen U, Sidhu H, Honow R, von Unruh G, Hesse A. The role of Oxalobacter formigenes colonization in calcium oxalate stone disease. *Kidney international*. Jun 2013;83(6):1144-1149.
253. Million M, Maraninchi M, Henry M, et al. Obesity-associated gut microbiota is enriched in Lactobacillus reuteri and depleted in Bifidobacterium animalis and Methanobrevibacter smithii. *International journal of obesity*. Jun 2012;36(6):817-825.
254. Mathur R, Amichai M, Chua KS, Mirocha J, Barlow GM, Pimentel M. Methane and hydrogen positivity on breath test is associated with greater body mass index and body fat. *The Journal of clinical endocrinology and metabolism*. Apr 2013;98(4):E698-702.
255. Pimentel M, Mayer AG, Park S, Chow EJ, Hasan A, Kong Y. Methane production during lactulose breath test is associated with gastrointestinal disease presentation. *Digestive diseases and sciences*. Jan 2003;48(1):86-92.
256. Pique JM, Pallares M, Cuso E, Vilar-Bonet J, Gassull MA. Methane production and colon cancer. *Gastroenterology*. Sep 1984;87(3):601-605.
257. Weaver GA, Krause JA, Miller TL, Wolin MJ. Incidence of methanogenic bacteria in a sigmoidoscopy population: an association of methanogenic bacteria and diverticulosis. *Gut*. Jun 1986;27(6):698-704.
258. Peled Y, Weinberg D, Hallak A, Gilat T. Factors affecting methane production in humans. *Gastrointestinal diseases and alterations of colonic flora. Digestive diseases and sciences*. Mar 1987;32(3):267-271.
259. Scanlan PD, Shanahan F, Marchesi JR. Human methanogen diversity and incidence in healthy and diseased colonic groups using mcrA gene analysis. *BMC Microbiol*. 2008;8:79.
260. Angelakis E, Armougom F, Million M, Raoult D. The relationship between gut microbiota and weight gain in humans. *Future microbiology*. Jan 2012;7(1):91-109.
261. Eckburg PB, Bik EM, Bernstein CN, et al. Diversity of the human intestinal microbial flora. *Science (New York, N.Y.)* Jun 10 2005;308(5728):1635-1638.
262. Moore WE, Moore LH. Intestinal floras of populations that have a high risk of colon cancer. *Applied and environmental microbiology*. Sep 1995;61(9):3202-3207.
263. Warren RL, Freeman DJ, Pleasance S, al e. Co-occurrence of anaerobic bacteria in colorectal carcinomas. *Microbiome* 2013;1:16.
264. McCoy AN, Araujo-Perez F, Azcarate-Peril A, Yeh JJ, Sandler RS, Keku TO. Fusobacterium is associated with colorectal adenomas. *PLoS One*. 2013;8(1):e53653.



265. Ahn J, Sinha R, Pei Z, et al. Human gut microbiome and risk for colorectal cancer. *J Natl Cancer Inst.* Dec 18 2013;105(24):1907-1911.
266. Swidsinski A, Dorffel Y, Loening-Baucke V, et al. Acute appendicitis is characterised by local invasion with *Fusobacterium nucleatum/necrophorum*. *Gut.* Jan 2011;60(1):34-40.
267. Andersen RN, Ganeshkumar N, Kolenbrander PE. *Helicobacter pylori* adheres selectively to *Fusobacterium* spp. *Oral microbiology and immunology.* Feb 1998;13(1):51-54.
268. Chen Y, Yang F, Lu H, et al. Characterization of fecal microbial communities in patients with liver cirrhosis. *Hepatology (Baltimore, Md.).* Aug 2011;54(2):562-572.
269. Strauss J, Kaplan GG, Beck PL, et al. Invasive potential of gut mucosa-derived *Fusobacterium nucleatum* positively correlates with IBD status of the host. *Inflammatory bowel diseases.* Sep 2011;17(9):1971-1978.
270. Allen-Vercoe E, Strauss J, Chadee K. *Fusobacterium nucleatum*: an emerging gut pathogen? *Gut microbes.* Sep 1 2011;2(5):294-298.
271. Derrien M, Vaughan EE, Plugge CM, de Vos WM. *Akkermansia muciniphila* gen. nov., sp. nov., a human intestinal mucin-degrading bacterium. *International journal of systematic and evolutionary microbiology.* Sep 2004;54(Pt 5):1469-1476.
272. Everard A, Belzer C, Geurts L, et al. Cross-talk between *Akkermansia muciniphila* and intestinal epithelium controls diet-induced obesity. *Proceedings of the National Academy of Sciences of the United States of America.* May 28 2013;110(22):9066-9071.
273. Papa E, Docktor M, Smillie C, et al. Non-invasive mapping of the gastrointestinal microbiota identifies children with inflammatory bowel disease. *PLoS One.* 2012;7(6):e39242.
274. Krznaric Z, Vranesic Bender D, Kunovic A, Kekez D, Stimac D. Gut microbiota and obesity. *Digestive diseases (Basel, Switzerland).* 2012;30(2):196-200.
275. Murphy EF, Cotter PD, Hogan A, et al. Divergent metabolic outcomes arising from targeted manipulation of the gut microbiota in diet-induced obesity. *Gut.* Feb 2013;62(2):220-226.
276. De Bandt JP, Waligora-Dupriet AJ, Butel MJ. Intestinal microbiota in inflammation and insulin resistance: relevance to humans. *Current opinion in clinical nutrition and metabolic care.* Jul 2011;14(4):334-340.
277. De Minicis S, Rychlicki C, Agostinelli L, et al. Dysbiosis contributes to fibrogenesis in the course of chronic liver injury in mice. *Hepatology.* May 2014;59(5):1738-1749.
278. Remely M, Aumueller E, Jahn D, Hippe B, Brath H, Haslberger AG. Microbiota and epigenetic regulation of inflammatory mediators in type 2 diabetes and obesity. *Benef Microbes.* Mar 2014;5(1):33-43.
279. Martinez I, Lattimer JM, Hubach KL, et al. Gut microbiome composition is linked to whole grain-induced immunological improvements. *Isme J.* Feb 2013;7(2):269-280.
280. Wolf KJ, Lorenz RG. Gut Microbiota and Obesity. *Current obesity reports.* Mar 1 2012;1(1):1-8.
281. Abubakar I, Irvine L, Aldus CF, et al. A systematic review of the clinical, public health and cost-effectiveness of rapid diagnostic tests for the detection and identification of bacterial intestinal pathogens in faeces and food. *Health Technol Assess.* Sep 2007;11(36):1-216.
282. Liu J, Gratz J, Maro A, et al. Simultaneous detection of six diarrhoea-causing bacterial pathogens with an in-house PCR-luminex assay. *Journal of clinical microbiology.* Jan 2012;50(1):98-103.
283. Wiemer D, Loderstaedt U, von Wulffen H, et al. Real-time multiplex PCR for simultaneous detection of *Campylobacter jejuni*, *Salmonella*, *Shigella* and *Yersinia* species in fecal samples. *Int J Med Microbiol.* Nov 2011;301(7):577-584.
284. Cunningham SA, Sloan LM, Nyre LM, Vetter EA, Mandrekar J, Patel R. Three-hour molecular detection of *Campylobacter*, *Salmonella*, *Yersinia*, and *Shigella* species in feces with accuracy as high as that of culture. *Journal of clinical microbiology.* Aug 2010;48(8):2929-2933.
285. Takeda T, Bhattacharya SK, Nair GB. Cooperative studies on diarrheal diseases. *Acta Paediatr Jpn.* Dec 1993;35(6):583-590.
286. Kim DH, Lee BK, Kim YD, Rhee SK, Kim YC. Detection of representative enteropathogenic bacteria, *Vibrio* spp., pathogenic *Escherichia coli*, *Salmonella* spp., *Shigella* spp., and *Yersinia enterocolitica*, using a virulence factor gene-based oligonucleotide microarray. *Journal of microbiology (Seoul, Korea).* Oct 2010;48(5):682-688.
287. Kallinowski F, Wassmer A, Hofmann MA, et al. Prevalence of enteropathogenic bacteria in surgically treated chronic inflammatory bowel disease. *Hepatogastroenterology.* Sep-Oct 1998;45(23):1552-1558.
288. Kerckhoffs AP, Ben-Amor K, Samsom M, et al. Molecular analysis of faecal and duodenal samples reveals significantly higher prevalence and numbers of *Pseudomonas aeruginosa* in irritable bowel syndrome. *Journal of medical microbiology.* Feb 2011;60(Pt 2):236-245.
289. Amin OM. Seasonal prevalence of intestinal parasites in the United States during 2000. *The American journal of tropical medicine and hygiene.* Jun 2002;66(6):799-803.
290. Lacy BE, Rosemore J. *Helicobacter pylori*: ulcers and more: the beginning of an era. *The Journal of nutrition.* Oct 2001;131(10):2789S-2793S.



Gut Microbiota Clinical Association Summary*

	IBS	Inflammation	Immune Modulation	Metabolic Disorders	Autism	Diversity Association
Bacteroidetes Phylum						
<i>Bacteroides-Prevotella</i> Group	L	L/H	L	L/H		LD
<i>B. vulgatus</i>	L		L/H	L	H	LD
<i>Barnesiella</i> spp.			L			
<i>Odoribacter</i> spp.		L				
<i>Prevotella</i> spp.		H	H	H	L	LD
Firmicutes phylum						
<i>Anaerotruncus colihominis</i>				L		HD
<i>Butyrivibrio crossotus</i>				L		HD
<i>Clostridium</i> spp.		L/H			H	LD
<i>Coprococcus eutactus</i>	L				L	HD
<i>Faecalibacterium prausnitzii</i>		L/H		L		HD
<i>Lactobacillus</i> spp.	L/H			H		HD
<i>Pseudoflavonifractor</i> spp.						HD
<i>Roseburia</i> spp.	L	L	L	L		HD
<i>Ruminococcus</i> spp.	L/H	H			H	LD
<i>Veillonella</i> spp.	L/H				L	HD
Actinobacteria phylum						
<i>Bifidobacterium</i> spp.	L	L	L	L/H	L	HD
<i>B. longum</i>				H		HD
<i>Collinsella aerofaciens</i>	L	H		H		
Proteobacteria phylum						
<i>Desulfovibrio piger</i>	H	H		L	H	
<i>Escherichia coli</i>	H	H		L/H		
<i>Oxalobacter formigenes</i>				L		
Euryarchaeota phylum						
<i>Methanobrevibacter smithii</i>	L/H	L/H		L/H		HD
Fusobacteria phylum						
<i>Fusobacterium</i> spp.		H	H	H		
Verrucomicrobia phylum						
<i>Akkermansia muciniphila</i>		L		L	L	HD

Key

L = Low

H = High

LD = Robust levels of this organism associated with Low Diversity of gut bacteria

HD = Robust levels of this organism associated with High Diversity of gut bacteria

*The literature-based clinical associations in this chart are not intended to indicate diagnostic patterns.

Asheville

63 Zillicoa Street
Asheville, NC 28801

800.522.4762
www.gdx.net

Atlanta

3425 Corporate Way
Duluth, GA 30096

800.522.4762
www.gdx.net

London

Parkgate House
356 West Barnes Lane
New Malden Surrey, KT3 6NB

020.8336.7750
www.gdxuk.net



GI *fx* GI Effects Stool Profiles®

GI Effects Stool Profiles

- See link below* for the following:
 - > # 2200 GI Effects
 - > # 2205 GI Effects Microbial Ecology Profile

* <http://www.gdx.net/tests/alphabetical#G>